TUMARK® Professional


REF 271557 271559 271560 271565

Tumark® Professional X-Shape

REF 271549 271551

Tumark® Professional Q-Shape

REF 271548 271555

GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MODE D'EMPLOI
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ISTRUZIONI PER L'USO
INSTRUÇÕES DE UTILIZAÇÃO
INSTRUCCIONES DE USO
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
NÁVOD K POUŽITÍ
UPUTE ZA UPORABU
ИНСТРУКЦИИ ЗА УПОТРЕБА
BRUGSANVISNING
KASUTUSJUHEND
KÄYTTÖOHJEET
HASZNÁLATI UTASÍTÁSOK
LIETOŠANAS INSTRUKCIJA
NAUDOJIMO INSTRUKCIJA
INSTRUCTIES VOOR GEBRUIK
BRUKSANVISNING
INSTRUKCJA UŻYCIA
INSTRUCTIUNI DE UTILIZARE
NÁVOD NA POUŽITIE
BRUKSANVISNING

INHALT / CONTENT

Contenu / Περιεχόμενο / Contenuto / Índice / Содержание / Obsah / Sadržaj / съдържание / indhold / Sisu / pitoisuus / tartalom
Satur / Turinys / Inhoud / innhold / Treść / conținut / obsah / innehåll:

| | |
|--------------------------|----|
| INHALT / CONTENT | 2 |
| DE - DEUTSCH | 3 |
| EN - ENGLISH | 6 |
| FR - FRANÇAIS | 9 |
| EL - ΕΛΛΗΝΙΚΑ | 12 |
| IT - ITALIANO | 15 |
| PT - PORTUGUÊS | 18 |
| ES - ESPAÑOL | 21 |
| RU - РУССКИЙ | 24 |
| CZ - ČESKY | 27 |
| HR - HRVATSKI | 30 |
| BG - БЪЛГАРСКИ | 33 |
| DA - DANSK | 36 |
| ET - EESTI | 39 |
| FI - SUOMI | 42 |
| HU - MAGYAR | 45 |
| LT - LIETUVIŲ K. | 48 |
| LV - LATVIEŠU | 51 |
| NL - NEDERLANDS | 54 |
| NO - NORSK | 57 |
| PL - POLSKI | 60 |
| RO - ROMÂNĂ | 63 |
| SK - SLOVENSKY | 66 |
| SV - SVENSKA | 69 |
| SYMBOLE / SYMBOLS | 72 |
| INFO | 78 |

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
0482

SOMATEX®
MEDICAL TECHNOLOGIES GMBH

EN - ENGLISH**Read instructions before use****Keep for future reference****Important Information:**

Read this instruction manual thoroughly and be familiar with its contents prior to use of the *Tumark® Professional*. Failure to read the entire manual and familiarize yourself with all instructions before using the *Tumark® Professional* is unsafe and can result in life-threatening or severe injury to the patient or user and to damage or malfunction of the device.

Intended use and indications for use:

The *Tumark® Professional* is intended for the percutaneous marking of soft tissue, such as breast tissue. Among its areas of application are marking of lesions prior to or during chemotherapy, marking the site of a biopsy or of a removed tumor, or prior to radiotherapy for better treatment planning.

Contraindications:

- The *Tumark® Professional* is not intended for use except as indicated above.
- The *Tumark® Professional* cannula is not suitable for use with magnetic resonance imaging (MRI).
- The *Tumark® Professional* is contraindicated in patients with severe nickel allergy.

A summary of safety and clinical performance in accordance with Article 32 of Regulation 2017/745 (implantable devices) is available on www.somatex.de.

Warnings:

- Only qualified physicians with knowledge, experience and training in percutaneous soft tissue marking shall use the *Tumark® Professional*.
- This manual does not include descriptions or instructions for surgical techniques. It is the responsibility of the physician performing any procedure to determine the appropriateness of the procedure to be performed and of the use of this device and to determine the specific technique for each patient.
- When implanting a clip marker near a breast implant, handle with care to avoid puncturing the breast implant.
- The *Tumark® Professional* should only be used if the indicator on the packaging is green, only before the expiry date and only if the packaging is unopened and undamaged. Product sterility can only be guaranteed, if these criteria are met. If the indicator is not green, if the expiry date is exceeded, or if the packaging is damaged or opened before use, the product should not be used and the distributor or manufacturer, SOMATEX, should be contacted.
- The product is intended for single use only: DO NOT reuse or resterilize.
- When using a positioning needle, the *Tumark® Professional* must be checked for compatibility in advance. The bevelled *Tumark® Professional* cannula tip opening should protrude fully out of the positioning needle, and the user should be able to gauge this protrusion in order to be able to apply the clip marker safely and not place it too far into the tissue.

Precautions:

- The *Tumark® Professional* clip marker is made from a nickel-titanium alloy (Nitinol), which is why the product is contraindicated in patients with severe nickel allergy.
- Make sure that the slide button remains in the retracted position while the cannula is being put in position.
- The clip marker must be placed by pushing the slide button forward as far as possible to the stop position.
- There is risk of injury due to the sharp cannula tip. Use care especially when unpacking the cannula.
- The cannula of *Tumark® Professional* is NOT made of MRI-compatible metals. NOT suitable for MRI safety area. Danger of injury!
- Pay attention to the dimensions of the clip marker in relation to the size of the tissue area being marked (see Product Description).
- In rare cases the expansion of the clip marker may be delayed. Visibility in radiological imaging might be compromised until full expansion.

Information about materials used:

The implantable clip marker is made from a nickel-titanium alloy (Nitinol).


Unlike the clip marker, the *Tumark® Professional* cannula application device is not suitable for use in an MRI scanner.

MR unsafe

MRI Safety Information clip marker:

The clip marker is conditionally MR safe. A patient can safely undergo an MRI procedure with the clip marker under the following conditions:

- static magnetic field up to 3.0 Tesla with
- theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg.

0482

SOMATEX®
MEDICAL TECHNOLOGIES GMBH

Non-clinical tests were performed on the following systems:

- 3 Tesla Siemens Magnetom Trio (Siemens Medical, Erlangen, Germany) MRI with software Numaris 4, syngo MR (Version "B14, N4_VB15A_LATEST_20070519_P12");
- 3 Tesla Siemens Magnetom Skyra (Siemens Medical, Erlangen, Germany) MRI with software Numaris 4, syngo MR (Version "D11, N4_VB11D_LATEST_20110129_P3");
- 1.5 Tesla Siemens Magnetom Avanto (Siemens Medical, Erlangen, Germany) MRI with software Numaris 4, syngo MR (Version "B13, N4_VB13A_LATEST_20060607_P29").

Under the scanning conditions defined above, it is expected that clip marker of *Tumark® Professional* will produce the following maximum RF-related temperature rise:

- at 1.5 Tesla: < 1,0 °C (2 W/kg SAR) after 20 min of continuous scanning,
- at 3.0 Tesla: < 1,0 °C (2 W/kg SAR) after 20 min of continuous scanning.


Under the scanning conditions defined above, it is expected that clip marker of *Tumark® Professional* will produce the following image artefacts:

- U-Shape at 1.5 Tesla: 7.1 mm spin echo sequence; 6.5 mm gradient echo sequence;
- U-Shape at 3.0 Tesla: 9.0 mm spin echo sequence; 6.5 mm gradient echo sequence;
- X-Shape at 1.5 Tesla: 4.7 mm spin echo sequence; 3.1 mm gradient echo sequence;
- X-Shape at 3.0 Tesla: 2.7 mm spin echo sequence; 3.7 mm gradient echo sequence;
- Q-Shape at 1.5 Tesla: 4.3 mm spin echo sequence; 3.2 mm gradient echo sequence;
- Q-Shape at 3.0 Tesla: 4.7 mm spin echo sequence; 4.0 mm gradient echo sequence.

Do not expose the implanted clip marker to unconventional and non-standardized MRI techniques other than the ones listed above, because it has NOT BEEN TESTED for that purpose.

Device Description:

This is a sterile product for single use only and consists of a non-absorbable nickel-titanium clip marker (1), an introducer cannula (2) and a plastic handle. When new and unopened, the clip marker is contained within the cannula. The cannula tip is bevelled to help insertion, has markings 1 cm apart for measuring the depth of penetration, and a textured surface behind the cannula tip. The handle is equipped with a slide button (3) which allows one handed placement of the marker by pressing it forward. A safety catch system prevents the slide button from inadvertently moving forward and therefore prevents premature deployment of the marker. The clip markers for the *Tumark® Professional* are available in three different shapes: X-shape (4), U-shape (1 and 5) and Q-shape (6). The symbol of the clip marker shape is depicted on the handle (7). For details please refer to order numbers.

Directions for Use:

1. Before opening, make sure that the packaging is not already open or damaged, that the indicator on the packaging is green, and that it is within the expiry date.
2. Disinfect the puncture area and cover the area around it with sterile drapes if required.
3. Use suitable imaging methods (ultrasound, mammography) to identify the target area. NOTE: the *Tumark® Professional* cannula is not suitable for the MRI safety zone.
4. Open the packaging and remove the product from packaging.
5. Remove the cannula protection hose from the outer cannula by twisting.
6. Use the cannula (2) to puncture the target area, and insert into the tissue. The depth of insertion can be read from the markings on the cannula when positioning the cannula tip.
7. Check the position of the cannula tip using suitable imaging techniques, and adjust if appropriate.
8. Place the clip marker (1) by pushing the slide button (3) forward as far as it will go.
9. Verify and record the position of the clip marker (1).
10. Remove the cannula (2).
11. Treat the wound.

0482

12. After use: dispose the application device properly, following internal guidelines if appropriate; however, at least one suitable container intended for contaminated cannulas should be provided to ensure safe disposal.

Warning:

The company SOMATEX does not assume any liability for the use of this product or its components in case of re-sterilisation or reuse. This product may not be reused after a single application. The quality of the materials, coats and adhesive joints could degrade. Safe use is not guaranteed any longer. The product that is already used once is not designed for the required cleaning and sterilization processes. The sterility of the reprocessed disposable products is therefore not guaranteed. The risk of unwanted injuries and infections, especially cross-infections between patient and medical staff inappropriately increases.

Storage Instructions:

Keep away from sunlight and heat.

Store in a dry area at a temperature of 5 – 30 °C.

Any serious events that occur in relation to the product should be reported to SOMATEX Medical Technologies GmbH as well as the competent national authority.

5°C 30°C

0482

SOMATEX®
MEDICAL TECHNOLOGIES GMBH

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