
Tumark[®] Vision

REF 271589 271590

GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MODE D'EMPLOI
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ISTRUZIONI PER L'USO
INSTRUÇÕES DE UTILIZAÇÃO
INSTRUCCIONES DE USO
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
NÁVOD K POUŽITÍ
UPUTE ZA UPORABU

INHALT / CONTENT

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EN - ENGLISH*Read instructions before use**Keep for future reference***Important Information:**

Read this instruction manual thoroughly and be familiar with its contents prior to use. Failure to read the entire manual and familiarize yourself with all instructions before using the Tumark® Vision is unsafe and can result in life-threatening or severe injury to the patient or user and to damage or malfunction of the device.

Indications:

The Tumark® Vision is intended for the percutaneous marking of soft tissue, particularly breast tissue. Among its areas of application are especially the marking of lesions prior to or during chemotherapy, pinpointing the site of a biopsy or the removal of a tumor, or prior to radiotherapy for better treatment planning.

Contraindications:

- The Tumark® Vision is not intended for use except as indicated above.
- The use of the Tumark® Vision may be contraindicated in patients who suffer from a severe nickel allergy.

Warnings:


- Only qualified physicians with knowledge, experience and training in percutaneous soft tissue marking shall use the Tumark® Vision.
- This manual does not include descriptions or instructions for surgical techniques. It is the responsibility of the physician performing any procedure to determine the appropriateness of the procedure to be performed and of the use of this device and to determine the specific technique for each patient.
- When implanting a marker near a breast implant, handle with care to avoid puncturing the implant
- The Tumark® Vision tissue site marking system is only sterile if the indicator on the package is green, if used before the expiration date and if package is unopened and undamaged. DO NOT use if indicator is not green, after the expiration date or if package is open or damaged.
- For single use only. DO NOT reuse or resterilize.

Precautions:

- Make sure that the slide button remains in the retracted position while the cannula is being put in position.
- The marker must be placed by pushing the slide button forward as far as possible to the stop position.
- Cannula tip is sharp. Use care especially when unpacking the cannula.
- The cannula of Tumark® Vision is NOT made of MRI-compatible metals. NOT suitable for MRI safety area. Danger of injury!
- Consider the size of the clip marker in relation to the area to be marked (see picture 1).
- In rare cases the expansion of the marker may be delayed. Visibility in radiological imaging might be compromised until full expansion.

Device Description:

The Tumark® Vision is a sterile, single use, preloaded tissue site marking system consisting of a non-absorbable nickel-titanium marker (1), an introducer cannula (2) and a plastic handle. The introducer cannula is designed with a beveled tip for introduction, 1 cm depth markings and a roughening before the cannula tip for better placement with ultrasound imaging. The handle is equipped with a slide button (3) which allows one handed placement of the marker by pressing it forward. A safety catch system prevents the slide button from inadvertently moving forward and therefore prevents premature deployment of the marker. The clip marker has a spherical shape. The clip shape is labelled on the handle.

Behaviour of the clip marker in Magnetic Resonance Tomography (MRT):

The Somatex Tumark® Vision as an application system for clip markers is not suitable for use in MRI.

The clip marker itself, which has already been placed inside a patient, can be exposed to a magnetic field of up to 3.0 Tesla, for example in follow-up examinations. With regard to interaction (displacement) there is no increased risks to a patient with a Tumark Vision® clip marker implant.

The behaviour of the clip marker in MRI was tested under the following conditions (not clinically):

- Static magnetic field of 1.5 Tesla and 3.0 Tesla
- Gradient field of 542 Gauß/cm max.

The temperature increase of the Tumark® Vision clip marker was assessed as *not clinically significant*.

The non-clinical testing was performed with 1.5 Tesla on a MRI of type *Magnetom Avanto (Siemens Medical, Erlangen, Germany)* with *Numaris 4, syngo MR B17* software. The clip marker image artefact amounted to 4.3 mm during the spin echo sequence and 5.2 mm during the gradient echo sequence.

The non-clinical testing was performed with 3 Tesla on a MRI of type *Magnetom Skyra (Siemens Medical, Erlangen, Germany)* with *Numaris 4, syngo MR D13* software. The clip maker image artefact amounted to 5.2 mm during the spin echo sequence and 5.8 mm during the gradient echo sequence.

Do not expose the implanted Tumark® Vision clip marker to unconventional and non-standardized MRI-techniques other than the ones listed above, because it has NOT BEEN TESTED for that purpose.

Directions for Use:

1. Prior to opening the pack, make sure that the indicator on the package is green and that the package has not been opened and/or damaged. In addition, check the sterilization expiration date.
2. Disinfect the puncture area and cover the area around it with sterile drapes if required.
3. Locate the target area by using appropriate imaging methods. **NOTE:** the Tumark® Vision cannula is not suitable for MRI scans.
4. Pull the protective tube from the cannula hub with a twisting movement.
5. Puncture the target area with the cannula (2) and insert the cannula into the tissue. The puncture depth can be measured by the marks on the cannula when the cannula tip is being positioned.
6. Confirm the needle placement with appropriate imaging systems. If necessary correct the placement.
7. Place the marker (1) by firmly pushing the slide-button (3) forward as far as possible to stop position.
8. Confirm and document the location of the marker (1).
9. Remove the cannula (2).
10. Treat the wound.
11. After the procedure, please ensure the appropriate disposal of the cannula (2) in the proper cannula container.

Warning:

The company SOMATEX does not assume any liability for the use of this product or its components in case of re-sterilisation or reuse. This product may not be reused after a single application. The quality of the materials, coats and adhesive joints could degrade. Safe use is not guaranteed any longer. The product that is already used once is not designed for the required cleaning and sterilisation processes. The sterility of the reprocessed disposable products is therefore not guaranteed. The risk of unwanted injuries and infections, especially cross-infections between patient and medical staff inappropriately increases.

Storage Instructions:

Keep dry.

Keep away from sunlight and heat (temperature between 5 – 30 °C).

This product is a medical device in accordance to directive 93/42/EEC.

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SOMATEX Medical Technologies GmbH