**MARROW CUT
BONE MARROW
SAFE CUT[®]**

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Important Information

Make yourself familiar with the entire content of these instructions for use and read them carefully before using the product.

Area of Use

The systems for bone marrow puncture Somatex[®] – Marrow Cut, Bone Marrow, and SAFE-CUT[®] – are used for the puncture and histological extraction of tissue from the bone marrow in the area of pelvic bone.

Contraindications

- Hemophilia or similar blood coagulation disorders
- Infection at the site of puncture

Warnings and Precautions

- The bone marrow puncture systems Somatex[®], described in these instructions for use, are not intended for sternal aspiration.
- For safe and effective application of bone marrow puncture systems, the physician must have relevant knowledge, experience and training in the use of this technique in patients.
- Data from clinical experiments indicates that it is absolutely necessary to follow Good Clinical Practice and preventive measures. Deep wound infections may occur as severe subsequent complications following application, and elimination thereof requires a major surgical procedure.
- These instructions are not intended as a description of techniques used during the intervention or surgical techniques. The performing physician is responsible for testing and considering the possibility of the safe application of these products and techniques for each individual patient.
- All components of the bone marrow puncture systems described herein remain sterile only until the first opening of and/or any damage to the packaging.
- The tools are intended explicitly for single use only. The product may not be re-used after a single application. The quality of materials, coatings and adhesive joints may deteriorate. A secure application is no longer guaranteed. The product is not designed for cleaning and sterilization processes after a single use. The sterility of reprocessed single-use is therefore not guaranteed. The risk of unintentional injuries and infections, especially cross-infection with a patient and medical personnel increases unmeasurably. The firm Somatex[®] – Medical Technologies – has no liability in the case of re-sterilization and/or repeated use of any tool.
- Caution, injury hazard! The individual components of these products have very sharp edges and can easily cause injury.
- Never use the product after the sterility expiration date indicated on the packaging.

Complications

- Additional bleeding and hematoma
- Infection
- Injury of the skin, soft tissues, nerves or adjacent organs
- Disorder of breathing functions during the use of analgesics or sedatives

Preventive Measures

- Bone marrow punctures should be performed in clinical settings only.
- Sterile disposable gloves should be worn.
- The area of puncture should be covered to create aseptic conditions.
- Appropriate measures should be taken prior to the puncture to assure biopsy sample for hemopathological evaluation.

Biopsy Application: (you should unconditionally follow the illustrations I and II on page 2)

1. Perform subcutaneous and periosteal local anesthesia at the puncture site.
2. When necessary, use a sharp scalpel to make an incision at the puncture site to facilitate penetration through the skin.
3. Prior to the beginning of the puncture, you can adjust the insertion depth stopper (3) on the cannule shaft as needed to ensure optical control of the penetration depth.
4. Insert the puncture needle (2) and advance it forward until you feel a noticeable resistance indicating that you reached the bone tissue.

-
5. Aim the needle and penetrate it into the crest of the pelvic bone under constant pressure using the turning motion (turning left/right).
 6. After reaching the bone marrow cavity (reduced resistance of penetration) pull the trocar out of the puncture needle by removing the closure lid (1).
 7. Optionally you can now attach a Luer-Lock syringe (10 ml) to the needle handle and aspirate blood from the bone marrow cavity.
 8. You can put the additional cap (5) on the needle handle to achieve better handling and to avoid further leakage of blood from the Luer-Lock connector.
 9. Now push the needle further into the bone marrow cavity by using rotating moves to cut the biopsy sample. The length of the biopsy sample can be checked by inserting the sliding biopsy sample gauge (7). Where necessary for this purpose, remove the additional cap (5) from the handle of the puncture needle (2) inserted into the crest of the pelvic bone and insert the sliding biopsy sample gauge (7) carefully into the cannule until you feel a noticeable resistance. Markers (8) on the body of the sliding biopsy sample gauge will indicate how long the biopsy sample is in the inner part of the puncture needle (the distance between the marker bands is 1 cm).
 10. Move the needle with a suitable strain back and forth, in order to separate the Biopsat from mark range and then remove the needle from the bones.
PRECAUTION: In this case avoid an unsuitable, excessive strain, in order to prevent the possible bending of the needle.
 11. Treat the wound.

Collection of the biopsy sample from the needle after the puncture

12. Remove the alternate cap (1) from the puncture needle (2).
13. Place the pusher guide (6) on the tip of the puncture needle in the direction of the arrow.
14. Insert the sliding biopsy sample gauge (7) into the conical orifice of the biopsy pusher guide and carefully remove the biopsy sample from the cannule in the direction from the distal to the proximal end.
15. Dispose of the used puncture needle in the appropriate manner.

When taking a biopsy sample using the puncture system of the Safe-CUT[®] type, you should proceed as follows:

(you should unconditionally follow illustration III on page 2)

- a) As soon as you have cut the biopsy sample and checked its sufficient length (see also step 9), remove the closing cap (5) and insert the LockWire (securing wire) (4) into the handle of the puncture cannule (2) inserted in the crest of the pelvic bone, until the handle of the LockWire (4) clearly snaps into its place in the handle of the puncture needle (2).
- b) After attaching the LockWire, rotate the entire puncture needle 2 times around the axis of the cannula by at least 360 degrees. When doing this, be careful not to change the position of the cannule.
- c) Then continue as described in point 10 of the previous section.

Storage instructions:

Store in a dark, dry place (humidity 30 – 65%) at a temperature between 41 and 86° F (5 – 30°C).

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