CONTENT

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INSTRUCTIONS FOR USE (English)

Important Information:

Read this instruction manual thoroughly and be familiar with its contents prior to use. Failure to read the entire manual and familiarize yourself with all instructions before using the Tumark® MRI is unsafe and can result in life threatening or severe injury to the patient or user and to damage or malfunction of the device.

Indications:

The Tumark® MRI is intended to attach a marker to soft tissue at the surgical site during an open or a percutaneous procedure. It is indicated for use to radiographically and radiologically mark the location in breasts tissue following an open or percutaneous procedure. It is intended to be used with magnetic resonance imaging (MRI) techniques.

Purpose of the Device:

The Tumark® MRI serves for marking of soft breast tissue which ensures visibility using magnetic resonance imaging (MRI) techniques and therefore lesion localization at a later date. The marker may be implanted to mark the location of a biopsy sampling point or the location of a removed tumor. In addition, it may be implanted in lesions prior to or during chemotherapy.

Contraindications:

- The Tumark® MRI is not intended for use except as indicated above.
- The use of the Tumark® MRI is contraindicated in patients who suffer from a severe nickel allergy.

Warnings:

- Only qualified physicians with knowledge, experience and training in percutaneous soft tissue marking should use the Tumark® MRI.
- Only qualified physicians with knowledge, experience and training in MR Imaging should use the Tumark® MRI in Magnetic Resonance Imaging (MRI) procedures.
- This manual does not include descriptions or instructions for surgical techniques. It is the responsibility of the physician performing any procedure to determine the appropriateness of the type of procedure to be performed and the use of this product and to determine the specific technique for each patient.
- DO NOT use the system in patients with breast implants.
- The Tumark® MRI tissue site marking system is only sterile if the indicator is green, if used before the expiration date and if package is unopened and undamaged. DO NOT use if indicator is not green, after the expiration date or if package is open or damaged.
- Single patient use only. DO NOT reuse or resterilize.
- Note that the size of the image artifact depends on the sequence and the orientation towards the BO field. Therefore the actual position of the cannula tip may differ from the position of the artifact tip.

Precautions:

- Make sure that the slide-button does not change its backward position while placing the cannula.
- The marker must be placed by pushing the slide-button forward as far as possible.
- Cannula tip is sharp. Use care especially while unpacking the cannula.
- **Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician.**
- The Tumark® MRI is MR Conditional. It can be safely scanned under the following conditions:
 - Static magnetic field of 1,5 – Tesla or less
 - Spatial gradient field of 260 Gauss / cm or less
 - Maximum averaged specific absorption rate (SAR) of 2 W/kg for 20 minutes of scanning.

In non-clinical testing, the Tumark® MRI produced a temperature rise of less than 1.42 °C at a maximum averaged specific absorption rate (SAR) of **2 W/kg**, as assessed by calorimetry for 20 minutes of MR scanning in a 1.5 Tesla Magnetom Avanto (Siemens Medical, Erlangen, Germany) MR scanner with software Numaris 4, syngo MR (version "B13, N4_VB13A_LATEST_20060607_P29").

The image artifact extends approximately 13.4 mm from the device when scanned in nonclinical testing using a spin echo sequence and approximately 8.9 mm from the device when scanned in nonclinical testing using a gradient echo sequence in a 1.5 Tesla Magnetom Avanto (Siemens Medical, Erlangen, Germany) MR system with software Numaris 4, syngo MR (version "B13, N4_VB13A_LATEST_20060607_P29").

STERILE EO

CE 0482

SOMATEX®
MEDICAL TECHNOLOGIES GMBH
MADE IN GERMANY

- Static magnetic field of 3 – Tesla or less
- Spatial gradient field of 380 Gauss / cm or less
- Maximum averaged specific absorption rate (SAR) of 2 W/kg for 20 minutes of scanning

In non-clinical testing the Tumark® MRI produced a temperature rise of less than 0.52 °C at a maximum whole body averaged specific absorption rate (SAR) of **2 W/kg**, as assessed by calorimetry for 20 minutes of MR scanning in a 3 Tesla Magnetom Trio (Siemens Medical, Erlangen, Germany) MR scanner with software Numaris 4, syngo MR, version is "B15, N4_VB15A_LATEST_20070519_P12".

The image artifact extends approximately 52.1 mm from the device when scanned in nonclinical testing using a spin echoe sequence and approximately 47,8 mm from the device when scanned in nonclinical testing using a gradient echoe sequence in a 3 Tesla Magnetom Trio (Siemens Medical, Erlangen, Germany) MR scanner with software Numaris 4, syngo MR, version is "B15, N4_VB15A_LATEST_20070519_P12".

- The Clip Marker for the Tumark® MRI System is MR Conditional. It can be safely scanned in the following conditions:
 - Static magnetic field of 1,5 – Tesla or less
 - Spatial gradient field of 260 Gauss / cm or less
 - Maximum averaged specific absorption rate (SAR) of 2 W/kg for 20 minutes of scanning.

In non-clinical testing, the Clip Marker for the Tumark® MRI System produced a temperature rise of less than 0.72 °C at a maximum averaged specific absorption rate (SAR) of **2 W/kg**, as assessed by calorimetry for 20 minutes of MR scanning in a 1.5 Tesla Magnetom Avanto (Siemens Medical, Erlangen, Germany) MR scanner with software Numaris 4, syngo MR (version "B13, N4_VB13A_LATEST_20060607_P29").

The image artifact extends approximately 7.1 mm from the device when scanned in nonclinical testing using a spin echoe sequence and approximately 6.5 mm from the device when scanned in nonclinical testing using a gradient echoe sequence in a 1.5 Tesla Magnetom Avanto (Siemens Medical, Erlangen, Germany) MR system with software Numaris 4, syngo MR (version "B13, N4_VB13A_LATEST_20060607_P29").


- Static magnetic field of 3 – Tesla or less
- Spatial gradient field of 260 Gauss / cm or less
- Maximum averaged specific absorption rate (SAR) of 2 W/kg for 20 minutes of scanning

In non-clinical testing the Clip Marker for the Tumark® MRI System produced a temperature rise of less than 0.61 °C at a maximum whole body averaged specific absorption rate (SAR) of **2 W/kg**, as assessed by calorimetry for 20 minutes of MR scanning in a 3 Tesla Magnetom Trio (Siemens Medical, Erlangen, Germany) MR scanner with software Numaris 4, syngo MR, version is "B15, N4_VB15A_LATEST_20070519_P12".

The image artifact extends approximately 9 mm from the device when scanned in nonclinical testing using a spin echoe sequence and approximately 6.5 mm from the device when scanned in nonclinical testing using a gradient echoe sequence in a 3 Tesla Magnetom Trio (Siemens Medical, Erlangen, Germany) MR scanner with software Numaris 4, syngo MR, version is "B15, N4_VB15A_LATEST_20070519_P12".

Device Description:

The Tumark® MRI is a sterile, single use, preloaded tissue site marking system consisting of a non-absorbable nickel-titanium marker (1), an introducer cannula (2) and a plastic handle. The introducer cannula is designed with a beveled tip for convenient introduction and 1 cm depth marks to aid in cannula placement. The handle is equipped with a slide-button (3) which allows for a one handed placement of the marker by pressing it forward. A safety catch system prevents the slide-button to inadvertently move forward and therefore prevents a premature deployment of the marker.

Directions for Use:

1. Prior to opening the package, make sure that the indicator on the package is green and that the package has not been opened and/or damaged. In addition, check the sterilization expiration date.
2. Open the package.
3. Disinfect the puncture area and cover it with sterile drapes.
4. Locate the target area by using appropriate MR imaging systems.
5. Puncture the target area with the cannula (2) and insert the cannula into the breast. The depth marks can assist you while positioning the cannula to the target.
6. Confirm the needle placement with appropriate MR imaging systems. If necessary correct the placement.
7. Place the marker (1) by firmly pushing the slide-button (3) forward as far as possible.
8. Confirm and document the location of the marker (1).
9. Remove the cannula (2).
10. Treat the wound.
11. After the procedure, please ensure the appropriate disposal of the cannula (2) in the proper cannula container.

Warning:

The company SOMATEX does not assume any liability for the use of this product or its components in case of re-sterilisation or reuse. This product may not be reused after a single application. The quality of the materials, coats and adhesive joints could degrade. Safe use is not guaranteed any longer. The product that is already used once is not designed for the required cleaning and sterilisation processes. The sterility of the reprocessed disposable products is therefore not guaranteed. The risk of unwanted injuries and infections, especially cross-infections between patient and medical staff inappropriately increases.

Storage Instructions:

Store in a dark, dry (humidity 30% - 65%) place with a temperature between 41° - 86° F (5° - 30°C).

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Somatex Medical Technologies GmbH