Embryo Transfer Set - S

REF 420322

GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MODE D'UTILISATION
ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ
ISTRUZIONI PER L'USO
INSTRUÇÕES DE UTILIZAÇÃO
MODO DE EMPLEO
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
NÁVOD K POUŽITÍ

INHALT / CONTENT

INHALT / CONTENT	2
DEUTSCH	3
ENGLISH.....	5
FRANÇAIS	7
ΕΛΛΗΝΙΚΑ	9
ITALIANO.....	11
PORTUGUÊS.....	13
ESPAÑOL	15
РУССКИЙ.....	17
ČESKY	19
SYMBOLE / SYMBOLS.....	21
INFO	24

STERILE R**CE**
0482**SOMATEX®**
MEDICAL TECHNOLOGIES GMBH

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ENGLISH**WARNING:**

Please read all the information carefully! Misuse can lead to serious complications.

Intended Purpose:

The SOMATEX® Embryo Transfer Set consists of a cervical catheter and an embryo transfer catheter, and serves the introduction of *in vitro* fertilised (IVF) embryos into the uterine cavity.

Product Description:

The SOMATEX® Embryo Transfer Set is a sterile disposable product (single-use) consisting of two low-friction coated catheters. The outer cervical catheter allows insertion through the cervix. The inner transfer catheter can be connected to a Luer Lock syringe via a handle.

Contraindications:

SOMATEX® catheters should not be used in patients with vaginal or intrauterine infection, sexually transmitted infection, current or recent pregnancy, recent uterine perforation, or in patients who have an IUD in the uterus.

Possible Known Complications:

During vaginal use, infections can occur through bacterial contamination (e.g. inflammatory pelvic diseases, urinary tract infections, cystitis and intrauterine infections). In addition, it is recommended that the catheter (as well as all other aids) is rinsed with sterile, embryo-compatible culture media and that the work is performed in a strictly sterile manner.

When inserting a catheter, bleeding may occur as a result of trauma. This has been associated with a lower pregnancy rate. Experience has shown that a simple and atraumatic transfer method offers the best chances of success.

Warnings and Precautions:

- The product may only be used by qualified personnel with appropriate knowledge, experience and training.
- These instructions are NOT intended to specify or suggest medical or surgical methods. The respective physician is responsible for determining the appropriate procedures and methods in connection with the use of this product.
- We recommend using only materials that have been tested for compatibility with mouse embryos.
Please note: The blastocyst rate which was determined for the SOMATEX® Embryo Transfer Set with single-cell mouse embryos is at least 80%. The LAL tests (endotoxin tests) were passed with a maximum of 20 EU (endotoxin units) per product. These tests are carried out in batches.
- The handling of biological samples and accessories must be carried out under observance of and in compliance with the relevant valid hygiene and safety regulations.
- HBV, HCV, HIV and other infectious diseases can be transmitted through contact with biological samples. In the case of exposure to biological samples and the resulting risk of infection, immediate medical attention is recommended.
- During removal or disposal, pay attention to the possible dripping of body fluids.
- SOMATEX® assisted reproduction products are expressly manufactured for single-use only. The respective product must be disposed of after use and not reused. The quality of the materials, coatings and adhesive bonds may deteriorate. Safe application is no longer guaranteed. The respective product is not designed for the necessary cleaning and sterilisation processes required following a single application. The sterility of the reprocessed disposable single-use products is therefore not guaranteed. The risk of inadvertent injury and infection, particularly cross-infection among patients and medical personnel, increases disproportionately. In the event of re-sterilisation or reuse, SOMATEX will decline any liability for the use of these products or individual product components.

Handling:

1. Before use, ensure that the package has not been opened and/or damaged. **WARNING:** If the packaging has already been opened or damaged, discard the product and choose one with the packaging intact.
2. Check the sterile date and open the sterile packaging.
3. Remove the product from the sterile packaging.
4. Position the patient in lithotomy or knee-chest position.
5. Introduce the vaginal speculum and expose the cervix.
6. Insert the outer cervical catheter.
7. Connect the transfer catheter to a syringe and aspirate the embryos into the distal catheter area.
8. Carefully insert the embryo-filled inner transfer catheter through the outer cervical catheter into the uterine cavity.
9. Then use the syringe to push the embryos into the uterine cavity.
10. Remove the catheters and check for any remaining embryos.
11. Dispose of the used SOMATEX® products in the disposal boxes provided according to the safety guidelines for your particular site.

Storage Instructions:

Store in a dry place.

Protect from direct sunlight and heat (temperature between 5–30 °C).

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SOMATEX Medical Technologies GmbH