
MRI Duo System

REF 601609 601611

GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MODE D'UTILISATION
ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ
ISTRUZIONI PER L'USO
INSTRUÇÕES DE UTILIZAÇÃO
MODO DE EMPLEO
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
NÁVOD K POUŽITÍ
UPUTE ZA UPORABU

INHALT / CONTENT

INHALT / CONTENT	2
ABBILDUNGEN/ PICTURES	3
DEUTSCH	4
ENGLISH	6
FRANÇAIS	8
ΕΛΛΗΝΙΚΑ	10
ITALIANO	12
PORTUGUÊS	14
ESPAÑOL	16
РУССКИЙ	18
ČESKY	19
HRVATSKI	21
SYMBOLE / SYMBOLS	24
INFO	27

STERILE EO**CE**
0482**SOMATEX®**
MEDICAL TECHNOLOGIES GMBH

ABBILDUNGEN/ PICTURES

Abbildung 1 / Picture 1:


Abbildung 2 / Picture 2:

Abbildung 3 / Picture 3:



Abbildung 4 / Picture 4:

DEUTSCH

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ENGLISH*Read carefully before use**Keep for future reference***Intended use:**

The product is intended as a preoperative marker of non-palpable suspected breast lesions to facilitate the intraoperative localization of the findings by the surgeon.

Product description:

The *MRI Duo System* consists of a hypodermic needle and a preloaded marking wire with distal double arches (Fig. 1).

Contraindications:

All contra-indications applicable to the relevant area of application, as known according to the rules of the art of medicine and anticipated for the use of cannulas and marker systems for preoperative marking of breast lesions, shall apply.

Possible known complications:

Wire dislocation, accidental cutting of wire, wire break, bleedings, infections, aesthetic complications

Warnings:

- Only qualified physicians with the required knowledge, experience and training shall use the product.
- Only qualified physicians with knowledge, experience and training in MR Imaging should use the *MRI Duo System* in Magnetic Resonance Imaging (MRI) procedures.
- This manual does not include descriptions or instructions for surgical techniques. It is the responsibility of the physician performing any procedure to determine the appropriateness of the procedure to be performed and of the use of this device and to determine the specific technique for each patient.
- The product is only sterile, if used before the expiration date and if package is unopened and undamaged. DO NOT use after the expiration date or if package is open or damaged.
- Note that the size of the image artifact depends on the sequence and the orientation towards the B₀ field. Therefore the actual position of the cannula tip may differ from the position of the artifact tip.
- Single patient use only. DO NOT reuse or resterilize.

Warning:

The company SOMATEX does not assume any liability for the use of this product or its components in case of re-sterilization or reuse. This product may not be reused after a single application. The quality of the materials, coats and adhesive joints can degrade. Safe use is not guaranteed any longer. The product that is already used once is not designed for the required cleaning and sterilization processes. The sterility of the reprocessed disposable products is therefore not guaranteed. The risk of unwanted injuries and infections, especially cross-infections between patient and medical staff inappropriately increases.

MRI Safety Information:

MRI Duo System is MR conditional. A patient with *MRI Duo System* can be safely scanned in an MR system meeting the following conditions:

- static magnetic field of **1.5 Tesla** and **3.0 Tesla** with
- theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR):
 - of 0.90 W/kg at 1.5 Tesla;
 - of 1.72 W/kg at 3.0 Tesla.

Non-clinical tests were performed on the following systems:

- 1.5 Tesla Philips Intera (Philips Medical Systems, Best, Netherlands) MRI with Software Release 12.6.1.4 2012-05-022;
- 3.0 Tesla Siemens Magnetom Trio (Siemens Medical, Erlangen, Germany) MRI with Software Numaris 4, syngo MR (Version "B17").

Under the scanning conditions defined above, it is expected that *MRI Duo System* will produce the following maximum RF-related temperature rise:

- at 1.5 Tesla: 6.0 °C (0.90 W/kg SAR) after 15 min of continuous scanning,
- at 3.0 Tesla: 6.0 °C (1.72 W/kg SAR) after 15 min of continuous scanning.

Under the scanning conditions defined above, it is expected that *MRI Duo System* will produce the following maximum image artefacts:

- at 1.5 Tesla: 8.9 mm spin echo sequence; 10.0 mm gradient echo sequence,
- at 3.0 Tesla: 52.1 mm spin echo sequence; 47.8 mm gradient echo sequence.

The following steps must be adhered to prior to the application:

1. Check that the package has not been opened or damaged before use; check the sterility date.
2. Open the package.

Note:

The shaft and the double arches of the marking wire are situated within the hypodermic needle in the delivered condition. Before using the marking wire, ensure that both of the first dark annular markings behind the needle handle are visible (Fig. 3, (3)). In this position, the arches of the wire are positioned completely within the needle.

Description of application:

3. Local disinfection and administration of local anesthetic in the insertion area.
4. Use the scalpel to carry out an incision in the puncture site to facilitate the insertion into the skin.
5. Needle insertion: insert the needle under guidance into the mamma such that the tip of the needle reaches the tumor itself or the tumor region.
6. After achieving the target point, you can place a rolled wire into the tumor area (Fig. 2, (2)). To achieve this, you should hold the needle position and introduce the wire into the cannula in the direction from the proximal to the distal point. Hold the wire about 1 cm away from the cannula hub and advance it into the cannula in several stages. Then the rolled wire should be completely deployed, as soon as the marking reaches the needle insertion (Fig. 2, (1)).
7. Prior to the removal of the needle, a proper positioning of the rolled wire can be inspected. Should the position be other than optimal, the wire can be withdrawn into the cannula. After rectifying the needle position, the wire can be deployed again according to the description in point 5. and 6. for preoperative marking.
8. Remove the needle from the mamma carefully.
9. Relax the compression of the breast, if necessary. Slide the opened fixation part over the wire almost to the surface of the skin of the breast (Important: The breast must be decompressed when this is done. There must still be some space between skin and fixation part). The wire must lie between the brackets (Fig. 4, (4)) of the fixation part. The wire is fixed through the closure being locked in place (Fig. 4, (5)). This measure prevents the risk of change in position of the wire along the distal end.
10. Treat the wound.
11. Dispose the product properly in a suitable container.

Storage instructions:

Protect from moisture. Keep dry.

Keep away from sunlight and heat (temperature 5 – 30 °C / 41 – 86 °F).

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SOMATEX Medical Technologies GmbH