
# ***Duo System Stereo***

**REF**

271991

271992

**GEBRAUCHSANWEISUNG**  
INSTRUCTIONS FOR USE  
MODE D'UTILISATION  
ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ  
ISTRUZIONI PER L'USO  
INSTRUÇÕES DE UTILIZAÇÃO  
MODO DE EMPLEO  
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ  
NÁVOD K POUŽITÍ  
UPUTE ZA UPORABU

## INHALT / CONTENT

INHALT / CONTENT .....	2
ABBILDUNGEN/ PICTURES:.....	2
DEUTSCH .....	3
ENGLISH.....	4
FRANÇAIS .....	5
ΕΛΛΗΝΙΚΑ .....	6
ITALIANO.....	7
PORTUGUÊS.....	8
ESPAÑOL .....	9
РУССКИЙ.....	10
ČESKY .....	11
HRVATSKI .....	12
SYMBOLE / SYMBOLS.....	13
INFO .....	16

## ABBILDUNGEN/ PICTURES:

Abbildung 1 / Figure 1:


Abbildung 2 / Figure 2:

Abbildung 3 / Figure 3:



Abbildung 4 / Figure 4:

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0482

**SOMATEX**<sup>®</sup>  
MEDICAL TECHNOLOGIES GMBH

**ENGLISH****Read carefully before use****Keep for future reference****Intended use:**

The product is intended as a preoperative marker of non-palpable suspected breast lesions to facilitate the intraoperative localization of the findings by the surgeon.

**Product description:**

The Duo System Stereo consists of a hypodermic needle and a preloaded marking wire with distal double arches (*figure 1*).

**Contraindications:**

All contra-indications applicable to the relevant area of application, as known according to the rules of the art of medicine and anticipated for the use of cannulas and marker systems for preoperative marking of breast lesions, shall apply.

**Possible known complications:**

Wire dislocation, accidental cutting of wire, wire break, bleedings, infections, aesthetic complications

**Warnings:**

- Only qualified physicians with the required knowledge, experience and training shall use the product.
- This manual does not include descriptions or instructions for surgical techniques. It is the responsibility of the physician performing any procedure to determine the appropriateness of the procedure to be performed and of the use of this device and to determine the specific technique for each patient.
- The product is only sterile, if used before the expiration date and if package is unopened and undamaged. DO NOT use after the expiration date or if package is open or damaged.
- Single patient use only. DO NOT reuse or resterilize.
- The product is NOT suitable for MRI use (magnetic resonance imaging)! Danger of injury!

**Note:**

The shaft and the double arches of the marking wire are situated within the hypodermic needle in the delivered condition. Before using the marking wire, ensure that both of the annular markings behind the cannula handle are visible (*figure 2, A*). In this position, the arches of the wire are positioned completely within the cannula.

**Caution:** Federal (USA) Law restricts this device to sale by or on the order of a physician.

**Application:**

1. Open the package.
2. Local disinfection and administration of local anesthetic in the insertion area.
3. Use the scalpel to carry out an incision in the puncture site to facilitate the insertion into the skin.
4. Needle insertion: insert the needle under sonographic or mammographic guidance into the mamma such that the tip of the needle reaches the tumor itself or the tumor region.
5. After achieving the target point, you can place the rolled wire into the tumor area (*figure 3*). To achieve this, you should hold the needle position and introduce the wire into the cannula in the direction from the proximal to the distal point. Hold the wire about 1 cm away from the cannula hub and advance it into the cannula in several stages. Then the rolled wire should be completely deployed, as soon as the marking reaches the cannula insertion (*figure 3, B*).
6. Prior to the removal of the cannula, a proper positioning of the rolled wire must be inspected. Should the position be other than optimal, the wire can be withdrawn into the cannula. After rectifying the cannula position, the wire can be deployed again - according to the description in point 5 - for preoperative marking.
7. Remove the cannula from the mamma carefully without dislocating the wire.
8. Relax the compression of the breast. Slide the opened fixation part, starting from the flat site (*figure 4, C*), over the wire almost to the surface of the skin of the breast.  
**Caution:** The breast must be decompressed, when this is done. There must still be some space between skin and fixation part.  
The fixation keel can be fixed to the wire by turning the screw cap. This measure prevents the risk of change in position of the wire along the distal end.
9. Treat the wound.
10. Dispose the product properly in a suitable container.

**Warning:**

The company SOMATEX does not assume any liability for the use of this product or its components in case of re-sterilization or reuse. This product may not be reused after a single application. The quality of the materials, coats and adhesive joints can degrade. Safe use is not guaranteed any longer. The product that is already used once is not designed for the required cleaning and sterilization processes. The sterility of the reprocessed disposable products is therefore not guaranteed. The risk of unwanted injuries and infections, especially cross-infections between patient and medical staff inappropriately increases.

**Storage:**

Protect from moisture. Keep dry.

Keep away from sunlight and heat (temperature 5 – 30 °C / 41 – 86 °F).

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**SOMATEX Medical Technologies GmbH**