# Aspi Cut


REF

181011

181013

181015

181016


**GEBRAUCHSANWEISUNG**  
**INSTRUCTIONS FOR USE**  
**MODE D'UTILISATION**  
**ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ**  
**ISTRUZIONI PER L'USO**  
**INSTRUÇÕES DE UTILIZAÇÃO**  
**MODO DE EMPLEO**  
**ИНСТРУКЦИЯ ПО**  
**ПРИМЕНЕНИЮ**  
**NÁVOD K POUŽITÍ**  
**UPUTE ZA UPORABU**

INHALT / CONTENT

INHALT / CONTENT .....	2
ABBILDUNGEN / PICTURES .....	2
DEUTSCH .....	3
ENGLISH .....	4
FRANÇAIS .....	5
ΕΛΛΗΝΙΚΑ .....	6
ITALIANO .....	7
PORTUGUÊS .....	8
ESPAÑOL .....	9
РУССКИЙ .....	10
ČESKY .....	11
HRVATSKI .....	12
SYMBOLE / SYMBOLS .....	13
INFO .....	16

ABBILDUNGEN / PICTURES

Photos / Εικόνες / Figure / Imagens / Imágenes / Иллюстрации / Obrázky / Slike

- 1 – **Kanülenrohr** / cannula tube / Tube de canule / Σωλήνας καθετήρα / tubo cannula / Tubo de cânula / Tubo de cánula / Канюля / Kanyla / cijev kanile
- 2 – **Mandrin** / mandrel / Stilet / Μανδρέλι / mandrino / Mandril / Mandril / Мандрен / Mandrén / stilet
- 3 – **Mandrinansatz** / mandrel hub / Moyeu du stilet / Υποδοχή μανδρελιού / perno mandrino / Base de mandril / Conector de mandril / Мундштук мандрена / Концовка мандрéну / priključak stileta
- 4 – **Ansatz** / hub / Moyeu / Υποδοχή / perno / Base / Conector / Мундштук / Концовка / priključak
- 5 – **Tiefenstopper** / Depth stopper / Butée de profondeur / Αναστολέας βύθισης / arresto di profondità / Limitador de profundidade / Торе profundo / Ограничитель глубины / Ηλιοβκόνά zarážka / graničnik za dubinu
- 6 – **Griff** / grip / Poignée / Λαβή / impugnatura / Pega / Mango / Рукоятка / Rukojeť / drška
- 7 – **Tiefenregulierung** / depth regulation / Régulation de profondeur / Ρύθμιση βάθους / regolazione profondità / Regulação de profundidade / Regulación profunda / Регулятор глубины / Regulace hloubky / regulator dubine

0482

**SOMATEX**<sup>®</sup>  
MEDICAL TECHNOLOGIES GMBH

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**ENGLISH***Read carefully before use**Keep until package is empty***INTENDED PURPOSE:**

Somatex *Aspi Cut* is used for cytological punctures in the sternal and/or posterior iliac crest area.

**CONTRAINDICATIONS:**

- Patients with a coagulation disorder under anticoagulant treatment require special medical assessment.
- Infections in the area of the puncture site

**POSSIBLE KNOWN COMPLICATIONS:**

- Subsequent bleeding and bruising
- Infections
- Damage to the skin, soft tissue, nerves or adjacent organs

**WARNINGS:**

- This product is for use by physicians only. Biopsy extraction requires explicit expertise in indications, contraindications, restrictions, typical findings and complications of a biopsy.
- For one-time use only. Do not resterilise. The risk of injuries caused by mechanical wear is increased considerably as is the risk of infections, in particular, cross-infections between patients and medical staff. In the event of resterilisation or reuse, SOMATEX will decline any responsibility for the use of this product or individual product components. This product must not be reused after a single application. The quality of the materials, coatings and adhesive joints may deteriorate. Safe application is no longer guaranteed. The product is not designed for the necessary cleaning and sterilisation processes after a single application. The sterility of the reprocessed disposable products is therefore not guaranteed.
- The contents are sterile. Do not use if the package is damaged. Please observe the expiration date.
- Do not use a damaged product.
- Handle with care. The cannula is sharp and there is a risk of injury.
- The product contains metal components. Not suitable for magnetic resonance tomography.

**WARNING:**

- The needle should be applied to the body under the control of an imaging process (sonography, tomography).
- Handle with care! Applying excessive force on the needle may bend it and interfere with the functioning. The needle is sharp and improper use can lead to injury.

**APPLICATION:**

1. Check the size of the cannula according to its colour coding. The colour coding table is in the appendix to these Instructions for Use.
2. Carefully remove the product from the sterile package. Carefully remove the protective sheath of the needle.
3. Perform local disinfection and local anaesthesia of the skin and periosteum in the area of the puncture site.
4. If necessary, make a stab incision over the puncture site with a sharp scalpel to facilitate penetration of the skin.
5. Use the depth regulation to configure the depth of penetration.
6. Advance the *Aspi Cut* puncture needle until resistance from bone tissue can be felt.
7. Insert the needle in a targeted manner, exerting continuous pressure whilst simultaneously rotating (clockwise / counter-clockwise) in the bone.
8. After reaching the medullary cavity (resistance to penetration ceases), remove the mandrel. First rotate the mandrel hub approximately 45° counter-clockwise.
9. Connect the needle with the Luer lock connection.
10. Aspirate the needle by pulling on the piston.
11. Disconnect the needle and secure the sample (possible additional sampling).
12. Remove the needle after adequate sampling by rotating and simultaneously pulling out from the bone.
13. Take care of the wound site.
14. Dispose of the used puncture needle.

**WARNING:** Handle with care. The cannula is sharp and there is a risk of injury.

**DISPOSAL:** The needle is sharp and can be contaminated. Please dispose of the product in a suitable container. Please observe the rules and directions for proper disposal.

**STORAGE:** Store in a dark, dry location at 5 – 30°C (41 - 86°F). Protect from heat and direct sunlight.

**ORDER:** Please use the Product Group Summary and Product Information on the last page for reorders.

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**SOMATEX Medical Technologies GmbH**