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ENGLISH**CAUTION:**

Read carefully before use! Misuse can lead to serious complications.

Intended use:

The SOMATEX® *Follicle Aspiration Set* is intended for the aspiration and irrigation of oocytes from the woman's ovaries.

Product description:

The SOMATEX® *Follicle Aspirations Set* consists of a single-lumen follicle aspiration needle as well as tube accessories for connection to a vacuum pump.

The needle is low-friction coated and echogenic.

Depending on the set variant, the tube accessory has a three-way stopcock for rinsing the follicle.

Contraindications:

Patients with vaginal or intrauterine infections, sexually transmitted infections or current pregnancy must not use this product. Use of the product is contraindicated in the case of an increased bleeding tendency or when taking anticoagulants.

Known possible complications:

- Pain and cramping in the lower abdomen may arise as a result of procedure. In the event of sudden severe pain, cramping and nausea, immediate diagnosis and treatment is indicated.
- Haematuria may arise due to penetration of the urine-filled bladder. This complication generally resolves itself after a period of 24 hours.
- When the cannula passes through the bladder, this may lead to an extravasation of urine in the abdominal cavity. Symptoms of this known complication must be monitored. In general, this is not connected either to pain or to undesired side effects.
- The procedure may lead to infections (e.g. pelvic inflammatory disease, urinary infections, cystitis and intrauterine infections, the formation of abscesses). It is recommended rinsing the cannula (as well as all other aids) with sterile, embryo-compatible culture media and to work in a strictly sterile manner.
- Ovarian torsion is a painful and rare, but possible, complication, which must be treated immediately in order to prevent a necrosis of the ovary.
- In rare cases, the formation of utero-vaginal fistulas has been described following the removal of the oocytes.
- The occurrence of vaginal bleeding in connection with transvaginal access by oocyte aspiration is known. The bleeding can usually be stopped quickly using compression. In rare cases, internal injuries to organs (e.g. perforation of the intestinal walls) and blood vessels with severe bleeding may occur. If these injuries remain unnoticed, this may lead to a threatening patient condition. It is imperative to ensure a corresponding follow-up monitoring of the patient.

Warnings and Precautions:

- The product may only be used by qualified specialists with the corresponding knowledge, experience and training.
- The product is designed to be connected to a vacuum pump. Aspiration using a syringe is not recommended.
- These instructions DO NOT serve to specify or suggest medicinal or surgical methods. The physician is responsible for determining the most appropriate procedures and methods in the application of this product.
- The physician is responsible for selecting the appropriate cannula diameter. SOMATEX® does not take responsibility for any injury to the patient in the event of selecting the wrong cannula.
- SOMATEX® recommends that only material which has been tested for compatibility with mouse embryos should be used.
Note: The blastocyst rate which has been determined with single-cell mouse embryos is at a minimum of 80%. The LAL tests (endotoxin tests) were passed with a maximum of 20 EU (endotoxin units). These tests are performed for each batch.
- This intervention should be performed under appropriate anaesthesia.
- There is a risk of injury due to sharp cannula tips. Handle with care.
- Biological samples and accessories must be handled in compliance with hygiene and safety.
- HBV, HCV, HIV and other infectious diseases may be transmitted through contact with biological samples. In the case of exposure to biological samples and a risk of infection, immediate precautionary medicinal measures are advised.
- Rinsing of the follicles may increase the duration of the intervention.
- This device has been manufactured for single use only. The product may not be re-used. The quality of materials, coatings and adhesive joints may deteriorate. Safe application is no longer guaranteed. The product is not designed to be cleaned and sterilised after single use. Therefore, the sterility of reprocessed single use products is not guaranteed. The risk of unintentional injuries and infections, especially cross-infection is unreasonably high in this case. In case of re-sterilisation or reuse, SOMATEX® does not accept any liability for the use of this product or product components.



Figure 2: Orientation of needle bevel:

Recommended Handling (Ultrasound-guided):

1. Before use, ensure that the packaging has not been opened and/or damaged.
2. Check the sterility date on the product.
CAUTION: If the packaging has already been opened or damaged or if the sterility date has expired, the product must not be used.
3. Position the patient in the lithotomy or knee-chest position.
4. When using a needle guide (not part of the set) do attach now to the ultrasound probe.
5. Insert the ultrasound probe vaginally and focus on the ovaries and oocytes.
6. Take the product from the sterile packaging.
7. Remove the protective tube.
8. Connect the tubing set to the vacuum pump and place the silicone plug on a suitable tube.
9. If necessary, the follicular aspiration needle can be fixed to the ultrasound probe holder.
10. Advance the needle, controlled by ultrasound, into the follicle.
 NOTE: The groove of the cannula adapter indicates the orientation of the needle edge (see Figure 2).
11. Aspirate or rinse the oocytes out of the ovaries under vacuum pressure (see figure 1 for the connection).
12. Repeat steps 10 & 11 to retrieve further oocytes. When focusing on the other ovary, first remove the needle and reposition the ultrasound probe.
13. Remove the product after the procedure.
14. Dispose the product correctly after use.

Storage Instructions:

Store in a dry place.
 Protect from direct sunlight and heat (temperature between 5 – 30 °C).

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SOMATEX Medical Technologies GmbH