
Q-Wire System

REF 271462 271463 271464 271465

GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MODE D'UTILISATION
ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ
ISTRUZIONI PER L'USO
INSTRUÇÕES DE UTILIZAÇÃO
MODO DE EMPLEO
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
NÁVOD K POUŽITÍ
UPUTE ZA UPORABU
HASZNÁLATI ÚTMUTATÓ
GEBRUIKSAANWIJZING
INSTRUKSJONER FOR BRUK
INSTRUKCJA UŻYCIA
INŠTRUKCIE NA POUŽÍVANIE
ANVÄNDNINGSSINSTRUKTIONER
BRUGSANVISNING

INHALT / CONTENT

CONTENU / ΠΕΡΙΕΧΟΜΕΝΟ / SODDISFARE / CONTENTE / CONTENIDO / СОДЕРЖАНИЕ / OBSAH / SADRŽAJ /
TARTALOM / INHOUD / INNHOLD / ZAWARTOŚĆ / INNEHÅLL / INDHOLD:

INHALT / CONTENT	2
DEUTSCH	3
ENGLISH	5
FRANÇAIS	7
ΕΛΛΗΝΙΚΑ	9
ITALIANO	11
PORTUGUÊS	13
РУССКИЙ	15
ESPAÑOL	17
ČESTINA	19
HRVATSKI	21
MAGYAR	23
NEDERLANDS	25
NORSK	27
POLSKI	29
SLOVENČINA	31
SVENSKA	33
DANSK	35
SYMBOLE / SYMBOLS	37
INFO	42

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Instructions for use

Important Information:

Read this instruction manual thoroughly and be familiar with its contents prior to use. Failure to read the entire manual and familiarize yourself with all instructions before using the **Q-Wire System** is unsafe and can result in life threatening or severe injury to the patient or user and to damage or malfunction of the device.

Indications:

The **Q-Wire System** is intended to mark non-palpable, suspicious soft tissue preoperatively to facilitate intraoperative localization of the findings by the surgeon. It is indicated for use for radiographically and ultrasound guided marking of the location in breast tissue before an open procedure. It is not indicated to be used with magnetic resonance imaging (MRI) techniques.

Purpose of the Device:

The **Q-Wire System** serves for marking of soft breast tissue which ensures radiographical and radiological visibility using ultrasound and mammography and therefore lesion localization. After achieving a target point with the puncture needle, the wire can be placed into a tumor area.

Contraindications:


- The **Q-Wire System** is not intended for use except as indicated above.

Warnings:

- Only qualified physicians with knowledge, experience and training in percutaneous soft tissue marking should use the **Q-Wire System**.
- This manual does not include descriptions or instructions for surgical techniques. It is the responsibility of the physician performing any procedure to determine the appropriateness of the type of procedure to be performed and the use of this product and to determine the specific technique for each patient.
- DO NOT use the system in patients with breast implants.
- The **Q-Wire System** tissue site marking system is only sterile if used before the expiration date and if package is unopened and undamaged. DO NOT use after the expiration date or if package is open or damaged.
- Single patient use only. DO NOT reuse or resterilize.
- The company SOMATEX does not assume any liability for the use of this product or its components in case of re-sterilization or reuse. This product may not be reused after a single application. The quality of the materials, coats and adhesive joints could degrade. Safe use is not guaranteed any longer. The product that is already used once is not designed for the required cleaning and sterilization processes. The sterility of the reprocessed disposable products is therefore not guaranteed. The risk of unwanted injuries and infections, especially cross-infections between patient and medical staff inappropriately increases.
- The product is NOT suitable for MRI use (magnetic resonance imaging)! Danger of injury!
- The product is NOT suitable in direct contact with electrosurgical instruments! Danger of spark discharging and injury!
- Care must be taken, when using electrosurgical instruments: the product can be damaged in direct contact. If damage to the marking wire occurs, a follow-up check and, if necessary, removal of wire fragments remaining in the body is required.

Precautions:

- Cannula tip is sharp. Use care especially while unpacking the cannula.
- Make sure that the marker wire is located inside of the cannula and that the curved part of the wire do not protrude out of the tip of the cannula. Both circular marks (1) of the wire must be visible. (see Picture 1)
- **Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician.**


Picture 1


Device Description:

The **Q-Wire System** is a sterile, single use, tissue site marking system consisting of a nickel-titanium marker wire (3) and a puncture needle (2) and 1 cm depth marks to aid in cannula placement (see Picture 2). The puncture needle is designed with a bevelled tip for convenient introduction. The marker wire is designed with two different markings, the first (1) marker allows the use to confirm that the curved part of the wire do not protrude out of the tip of the cannula prior to deployment. To avoid the risk of a change of the wire position in the distal direction a wire fixation (4) is usable.

Picture 2

Picture 3

Directions for Use:

1. Prior to opening the package, make sure that the package has not been opened and/or damaged. In addition, check the sterilization expiration date.
2. Open the package.
3. Disinfect the puncture area and administer local anesthetic in the insertion area.
4. Use the scalpel to carry out an incision in the puncture site to facilitate the insertion into the skin.
5. Before starting the puncture make sure that the curved part of the wire do not protrude out of the tip of the cannula. Both circular marks (1) of the wire must be visible.
6. Needle insertion: insert the needle under sonographic or mammographic guidance into the mamma such that the tip of the needle reaches the tumor itself or the tumor region.
7. After positioning the needle in the target area, the marker wire can be placed in the tumor area by holding the needle position and pushing the wire forward into the cannula. The marking wire tip is fully extended and in a curved shape once the second mark (1) on the marking wire gets invisible. (Picture 3)

Note: The wire should be pushed into the cannula only until the last visible marking on the wire (1) has been inserted into the cannula attachment!

IMPORTANT: Do not advance wire into the cannula any further past the second marking on the wire!

8. Prior to the removal of the needle, a proper positioning of the rolled wire can be inspected. Should the position be other than optimal, the wire can be withdrawn into the cannula. After rectifying the needle position, the wire can be deployed again according to the description in points 5 and 6 for preoperative marking.
9. Carefully remove the needle (2) from the mamma. If possible, bring the patient in an upright position. Relax the compression of the breast, if necessary. Slide the opened fixation part over the wire almost to the surface of the skin of the breast (**Caution:** The breast must be decompressed when this is done. There must still be some space between skin and fixation part). The wire must lie between the brackets (4a) of the fixation part. The wire is fixed through the closure being locked in place (4b). This measure prevents the risk of change in position of the wire along the distal end.
10. Treat the wound.
11. After the procedure, please ensure the appropriate disposal of the cannula (2) in the proper cannula container.

Storage Instructions:

Store in a dark, dry (humidity 30% - 65%) place with a temperature between 5° – 30°C (41° - 86° F).

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SOMATEX Medical Technologies GmbH