
Tuloc

REF 271660 271661 271663

GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MODE D'UTILISATION
ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ
ISTRUZIONI PER L'USO
INSTRUÇÕES DE UTILIZAÇÃO
MODO DE EMPLEO
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
NÁVOD K POUŽITÍ
UPUTE ZA UPORABU
KASUTUSJUHEND
KÄYTTÖOHJEET
HASZNÁLATI ÚTMUTATÓ
NAUDOJIMO INSTRUKCIJOS
LIETOŠANAS INSTRUKCIJA
GEBRUIKSAANWIJZING
INSTRUKSJONER FOR BRUK
INSTRUKCJA UŻYCIA
INŠTRUKCIE NA POUŽÍVANIE
ANVÄNDNINGSSINSTRUKTIONER
BRUGSANVISNING

INHALT / CONTENT

CONTENU / ΠΕΡΙΕΧΟΜΕΝΟ / SODDISFARE / CONTENTE / CONTENIDO / СОДЕРЖАНИЕ / OBSAH / SADRŽAJ / SISU / SISÄLTÖ / TARTALOM / TURINYS / SATURU / INHOUD / INNHOLD / ZAWARTOŚĆ / OBSAH / INNEHÅLL
INDHOLD:

INHALT / CONTENT	2
ABBILDUNGEN/ PICTURES	3
DEUTSCH	4
ENGLISH	5
FRANÇAIS	6
ΕΛΛΗΝΙΚΑ	7
ITALIANO	8
PORTUGUÊS	9
ESPAÑOL	10
РУССКИЙ	11
ČEŠTINA	12
HRVATSKI	13
EESTI	14
SUOMI	15
MAGYAR	16
LATVISKI	17
LIETUVIŲ	18
NEDERLANDS	19
NORSK	20
POLSKI	21
SLOVENČINA	22
SVENSKA	23
DANSK	24
SYMBOLE / SYMBOLS	25
INFO	30

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ENGLISH**Read carefully before use****Keep for future reference****Intended use:**

The product is intended as a preoperative marker of non-palpable suspected breast lesions to facilitate the intraoperative localization of the findings by the surgeon.

Product description:

The *TULOC* consists of a hypodermic needle and a preloaded marking wire with distal double arches (*Figure 1*). The marker wire has a marking (*Figure 1*) which provides the user with an orientation guide during the operation.

Contraindications:

All contra-indications applicable to the relevant area of application, as known according to the rules of the art of medicine and anticipated for the use of cannulas and marker systems for preoperative marking of breast lesions, shall apply.

Possible known complications:

Wire dislocation, accidental cutting of wire, wire break, bleedings, infections, aesthetic complications

Warnings:

- Only qualified physicians with the required knowledge, experience and training shall use the product.
- This manual does not include descriptions or instructions for surgical techniques. It is the responsibility of the physician performing any procedure to determine the appropriateness of the procedure to be performed and of the use of this device and to determine the specific technique for each patient.
- The product is only sterile, if used before the expiration date and if package is unopened and undamaged. DO NOT use after the expiration date or if package is open or damaged.
- The product is NOT suitable for MRI use (magnetic resonance imaging)! Danger of injury!
- Single patient use only. DO NOT reuse or resterilize.
- The product is NOT suitable in direct contact with electrosurgical instruments! Danger of spark discharging and injury!
- Care must be taken, when using electrosurgical instruments: the product can be damaged in direct contact. If damage to the marking wire occurs, a follow-up check and, if necessary, removal of wire fragments remaining in the body is required.

Warning:

The company SOMATEX does not assume any liability for the use of this product or its components in case of re-sterilization or reuse. This product may not be reused after a single application. The quality of the materials, coats and adhesive joints can degrade. Safe use is not guaranteed any longer. The product that is already used once is not designed for the required cleaning and sterilization processes. The sterility of the reprocessed disposable products is therefore not guaranteed. The risk of unwanted injuries and infections, especially cross-infections between patient and medical staff inappropriately increases.

The following steps must be adhered to prior to the application:

1. Check that the package has not been opened or damaged before use; check the sterility date.
2. Open the package.

Note:

The shaft and the double arches of the marking wire are situated within the hypodermic needle in the delivered condition. Before using the marking wire, ensure that both of the annular markings behind the cannula handle are visible (*Figure 3, C*). In this position, the arches of the wire are positioned completely within the cannula.

Application:

3. Local disinfection and administration of local anesthetic in the insertion area.
4. Use the scalpel to carry out an incision in the puncture site to facilitate the insertion into the skin.
5. Needle insertion insert the needle under sonographic or mammographic guidance into the mamma such that the tip of the needle reaches the tumor itself or the tumor region.
6. After achieving the target point, you can place a rolled wire into the tumor area (*Figure 2, B*). To achieve this, you should hold the needle position and introduce the wire into the cannula in the direction from the proximal to the distal point. Then the rolled wire should be completely deployed, as soon as the marking reaches the needle insertion (*Figure 2, A*). Note: The wire may be pushed into the cannula only until the last visible marking on the wire (*A*) has been inserted into the cannula attachment!
7. Prior to the removal of the needle, a proper positioning of the rolled wire can be inspected. Should the position be other than optimal, the wire can be withdrawn into the cannula. After rectifying the needle position, the wire can be deployed again according to the description in point 6 for preoperative marking.
8. Remove the needle from the mamma carefully.
9. Slide the opened fixation part over the wire almost to the surface of the skin of the breast. **Important:** The breast must be decompressed when this is done. There must still be some space between skin and fixation part. The wire must lie between the brackets (*Figure 4, D*) of the fixation part. The wire is fixed through the closure being locked in place (*Figure 4, E*). This measure prevents the risk of change in position of the wire along the distal end.
10. Treat the wound.
11. Dispose the product properly in a suitable container.

Storage:

Protect from moisture. Keep dry.

Keep away from sunlight and heat (temperature 5 – 30 °C / 41 – 86 °F).

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SOMATEX Medical Technologies GmbH