Duo System Premium

REF 271770 271771 271773

| CONTENT | Page |
|---|-------------|
| GEBRAUCHSANWEISUNG (GERMAN) _____ | 3 |
| INSTRUCTIONS FOR USE (ENGLISH) _____ | 4 |
| MANUEL D'UTILISATION (FRENCH) _____ | 5 |
| ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ (ΕΛΛΗΝΙΚΑ) _____ | 6 |
| ISTRUZIONI PER L'USO (ITALIAN) _____ | 7 |
| INSTRUÇÕES DE UTILIZAÇÃO (PORTUGUESE) _____ | 8 |
| РУКОВОДСТВО ПО ПРИМЕНЕНИЮ (РУССКИЙ) _____ | 9 |
| INDICACIONES DE USO (SPANISH) _____ | 10 |
| NÁVOD K POUŽITÍ (CZECH) _____ | 11 |
| UPUTE ZA UPORABU (CROATION) _____ | 12 |
| KASUTUSJUHEND (ESTONIAN) _____ | 13 |
| KÄYTTÖOHJEET (FINISH) _____ | 14 |
| HASZNÁLATI ÚTMUTATÓ (HUNGARIAN) _____ | 15 |
| NAUDOJIMO INSTRUKCIJOS (LATVIAN) _____ | 16 |
| LIETOŠANAS INSTRUKCIJA (LITHUANIAN) _____ | 17 |
| GEBRUIKSAANWIJZING (DUTCH) _____ | 18 |
| INSTRUKSJONER FOR BRUK (NORWEGIAN) _____ | 19 |
| INSTRUKCJA UŻYCIA (POLISH) _____ | 20 |
| INŠTRUKCIE NA POUŽÍVANIE (SLOVAK) _____ | 21 |
| ANVÄNDNINGSSINSTRUKTIONER (SWEDISH) _____ | 22 |
| BRUGSANVISNING (DANISH) _____ | 23 |
| SYMBOLE / SYMBOLS _____ | 24 |
| BESTELL INFORMATIONEN / REORDER INFORMATION _____ | 29 |

Abbildung 1 / Picture 1:



Abbildung 2 / Picture 2:

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0482

INSTRUCTIONS FOR USE (English)**Duo System Premium****Correctable Localization Kit****Intended use:**

The product is intended as a preoperative marker of non-palpable suspected breast lesions to facilitate the intraoperative localization of the findings by the surgeon.

Contraindications:

All contra-indications applicable to the relevant area of application, as known according to the rules of the art of medicine and anticipated for the use of cannulas and marker systems for preoperative marking of breast lesions, shall apply.

Possible known complications:

Wire dislocation, accidental cutting of wire, wire break, bleedings, infections, aesthetic complications

Warnings:

- Only qualified physicians with the required knowledge, experience and training shall use the product.
- This manual does not include descriptions or instructions for surgical techniques. It is the responsibility of the physician performing any procedure to determine the appropriateness of the procedure to be performed and of the use of this device and to determine the specific technique for each patient.
- The product is only sterile, if used before the expiration date and if package is unopened and undamaged. DO NOT use after the expiration date or if package is open or damaged.
- The product is NOT suitable for MRI use (magnetic resonance imaging)! Danger of injury!
- Single patient use only. DO NOT reuse or resterilize.
- The product is NOT suitable in direct contact with electrosurgical instruments! Danger of spark discharging and injury!
- Care must be taken, when using electrosurgical instruments: the product can be damaged in direct contact. If damage to the marking wire occurs, a follow-up check and, if necessary, removal of wire fragments remaining in the body is required.

Warning:

The company SOMATEX does not assume any liability for the use of this product or its components in case of re-sterilization or reuse. This product may not be reused after a single application. The quality of the materials, coats and adhesive joints can degrade. Safe use is not guaranteed any longer. The product that is already used once is not designed for the required cleaning and sterilization processes. The sterility of the reprocessed disposable products is therefore not guaranteed. The risk of unwanted injuries and infections, especially cross-infections between patient and medical staff inappropriately increases.

Before use:

1. Make sure that the packaging has not been opened and/or damaged; check sterilisation date.
2. Open the packaging.
3. When the wire arches are visible, insert into the cannula by pulling the lever (Fig. 1, 2) back to the stop.
4. Check that the arches (Fig. 1, 1) have completely disappeared into the cannula.

Application:

5. Disinfect locally and apply local anaesthetic around the entry point.
6. If necessary, use a scalpel to make a stab incision at the puncture point, to facilitate penetration of the skin.
7. Place tip of cannula in the area to be marked.
8. Unfold the arches by moving the lever (Fig. 1, 2) sideways to the left, then pushing it all the way forward (Fig. 1, 3).
9. Note position of the wire. If the arches are not in the correct position, retract them again by pulling back the lever (Fig. 1, 3). Move the cannula tip into the desired position, re-extend the arches and note the new position of the arches.
10. Before removing the cannula, turn lever (Fig. 1, 4) to the right until it clicks into place. The lever will not now retract the wire!
11. Carefully remove the cannula from the patient.
12. Relax the compression of the breast, if necessary. Slide the opened fixation part over the wire almost to the surface of the skin of the breast (important: The breast must be decompressed when this is done. There must still be some space between skin and fixation part). The wire must lie between the brackets (Fig. 2, 5a) of the fixation part. The wire is fixed through the closure being locked in place (Fig. 2, 5b). This measure prevents the risk of change in position of the wire along the distal end.
13. Dress the wound.
14. Dispose the product properly in a suitable container.

Storage:

Protect from moisture. Keep dry.

Keep away from sunlight and heat (temperature 5 – 30 °C / 41 – 86 °F).

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SOMATEX Medical Technologies GmbH