

# **Intra-Uterine-Shunt**

**REF**

410060

**GEBRAUCHSANWEISUNG**  
INSTRUCTIONS FOR USE  
INSTRUCCIONES DE USO  
ISTRUZIONI PER L'USO  
MODE D'UTILISATION  
INSTRUÇÕES DE UTILIZAÇÃO  
INSTRUKCJA UŻYCIA  
BRUGSANVISNING  
ANVÄNDNINGSINSTRUKTIONER  
INSTRUKSJONER FOR BRUK  
ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ  
GEBRUIKSAANWIJZING  
NÁVOD K POUŽITÍ

**INHALT / CONTENT**

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
0482

**SOMATEX®**  
MEDICAL TECHNOLOGIES GMBH

**EN - ENGLISH****Read instructions before use****Keep for future reference****Important Information:**

Read this instruction manual thoroughly and be familiar with its contents prior to use. Failure to read the entire manual and to familiarize yourself with all instructions before using the *intrauterine shunt* is unsafe and can result in life-threatening or severe injury to the patient or user and in damage or malfunction of the device.

**Intended use and indications for use:**

The *intrauterine shunt* is used for decompression of the bladder in fetuses diagnosed with narrowing of the lower urinary tract.

**Contraindications:**

- The *intrauterine shunt* is contraindicated for decompressing fluid accumulations outside the fetal bladder such as hydrothorax, pulmonary sequester, congenital pulmonary malformation of the airways (CPAM), diaphragmatic hernia, urinoma, ascites and cysts.
- The *intrauterine shunt* is contraindicated in patients with nickel allergy.

**Known potential complications:**

- Dislocation or obstruction of the shunt might occur. Therefore, regular checks to ensure its drainage capacity are recommended and multiple punctures for repeated shunt placement may be necessary.
- A dislocation of the shunt may require a surgical procedure to remove the shunt.
- The growing of the fetus may decrease the drainage function of the shunt.
- The following complications have been reported in connection with the procedure: chorioamnionitis, premature rupture of membranes, abortion, premature birth, need for induced premature birth and fetal death. Known maternal risks are injuries of the uterus with temporary infertility and complications related physical and psychological distress.

Fetal risks are perforations of inner organs, gastroschisis, abdominal wall hernia, ascitis and postnatal bladder dysfunction.

Like for every invasive intervention the risks of infections and related complications exist.

**A summary of safety and clinical performance in accordance with Article 32 of Regulation 2017/745 (implantable devices) is available on [www.somatex.de](http://www.somatex.de).**

**Warnings:**

- Only qualified physicians with knowledge, experience and training for the above mentioned indication, for prenatal therapy as well as for ultrasound-guided punctures (e.g. DEGUM II/III) shall use the *intrauterine shunt*. SOMATEX Medical Technologies GmbH company denies any responsibility in case of use for other indications.
- This manual does not include descriptions or instructions for surgical techniques. It is the responsibility of the physician performing any procedure to determine the appropriateness of the procedure to be performed and of the use of this device and to determine the specific technique for each patient.
- The decision about performing the procedure is the sole responsibility of the treating physician after informed consent of the patient.
- The postnatal removal of the shunt can be associated with surgery.
- The *intrauterine shunt* should only be used before the expiration date, and only if the packaging is unopened and undamaged. Product sterility can only be guaranteed if these criteria are met. If the expiration date is exceeded or if the packaging is damaged or opened before use, the product should not be used and the distributor or manufacturer, SOMATEX, should be contacted.
- The product is intended for single use only: DO NOT reuse or resterilize. The shunt may no longer be used, if there is a risk that it has been contaminated.
- Ensure a sterile environment and aseptic way of working.
- Before surgery, consider the size of the fetus and the space in the amniotic sac for drainage. Ensure that the length of the shunt (25 mm) is sufficient for the envisioned placement position.
- There is a risk that the proximal end of the shunt is accidentally placed in the wall of the uterus and that another shunt may need to be placed.

**Precautions:**

- The *intrauterine shunt* is partially made from a nickel-titanium alloy (Nitinol), which is why the product is contraindicated in patients with nickel allergy.
- There is risk of injury due to the sharp cannula tip. Use care especially when unpacking the cannula.
- The system does not allow an acquisition of samples of the drained fluid.
- The shunt has two X-ray markers located underneath the shunt brims. These are used for postnatal localization under X-ray imaging, if the shunt had been dislocated in the uterus.
- **Warning:** X-ray examinations during pregnancy are contraindicated.

5°C → 30°C

0482

**SOMATEX®**  
MEDICAL TECHNOLOGIES GMBH

- The cannula of the *intrauterine shunt* is NOT made of MRI-compatible metals. NOT suitable for MRI safety area. Danger of injury!
- The *intrauterine shunt* has NOT been tested for MR compatibility and is therefore not suitable for the MR security area. There could be a risk of injury in an MR application!

#### Information about materials used:

- The implantable shunt of the IUS system is made from a nickel-titanium alloy (Nitinol), which has been expanded with X-ray markers and a silicone membrane.

#### MRI Safety Information:


The *intrauterine shunt* is **not** suitable for use in an MRI scanner.

#### Product Description:

The system consists of the following components:

- |                             |                             |
|-----------------------------|-----------------------------|
| 1 – Protective tube shunt   | 7 – Ejector                 |
| 2 – Cannula with handle     | c – First marking           |
| 3 – Protective tube cannula | b – Second marking          |
| 4 – Preloader               | 8 – Protective tube ejector |
| 5 – Cap preloader (white)   | 9 – Shunt brim              |
| 6 – Shunt                   | 10 – Grip tab               |

#### Directions for Use:

##### A. Unpacking

- Before use, ensure that the package has not been opened and/or damaged and that it is within the expiration date.
- Open the packaging and remove the product from packaging.

##### B. Preloading of the shunt

Only directly before use, the shunt must be pulled into the cannula with the preloader.


The preloader (4) and the shunt (6) are connected with a blue thread in order to preload the shunt in the cannula (2).

0482


**SOMATEX®**  
MEDICAL TECHNOLOGIES GMBH

3. Hold the handle of the cannula (2) and pull on the handle of the preloader (4) with the white cap (5) (see Figure 1). **Caution:** Make sure that the protective tube (3) remains on the cannula (2). Do not loosen the white cap (5).


**Figure 1: Pull the handle of the preloader**

4. Pull the handle of the preloader until the black 10 mm marking (a) is visible behind the protective tube (3) (see Figure 2). **Caution:** Preload the shunt carefully and evenly to avoid damaging the shunt.


**Figure 2: Pull the handle of the preloader until the black marking (a) is visible.**

5. Loosen and remove the white cap (5) (see Figure 3). Below the cap, two free thread ends are visible.


**Figure 3: Loosen the cap**

6. Pull out the preloader (4) slowly and evenly (see Figure 4) and remove it so that both ends of the thread are completely free (see Figure 5). **Caution:** The free ends of the thread should not be fixed by hand!


**Figure 4: Remove preloader**

7. Cut off one of the two threads close to the protective tube (3) (i.e. 1 cm behind the needle tip) with a scalpel or scissors. **Caution:** Use a pair of sterile scissors or a sterile scalpel.


**Figure 5: Cut off one of the threads**

8. Pull the remaining thread until the thread is completely removed. **Caution:** The thread must be completely removed from the cannula. If this is not possible, because e.g. accidentally both ends of the thread had been cut, the system must be discarded.


*Figure 6: Pull out the remaining thread*

9. Remove protective tube of the shunt (1) from the proximal end of the cannula (see *Figure 6*). **Caution:** Remove the protective tube only when the shunt had been completely preloaded into the cannula. If needed, the grip tab (10) can be removed, too.
10. Proceed with the ejector and remove the protective tube (8) from the ejector (7) (see *Figure 7*). **Note:** Do not remove the protective tube (3) of the cannula (2) yet (see *Figure 6*).

*Figure 7: Ejector (7) with a protective tube (8)*

11. Insert the ejector (7) through the cannula handle (2) and advance until the first marking (c) is in front of the cannula handle (see *Figure 8*). **Caution:** Do not insert the preloader (4) into the cannula handle as this might lead to shunt damage.
- Note:** You should notice that approx. 1.5-2 cm before the marking (c), it becomes harder to push the ejector since the shunt is being pushed by the ejector in the cannula.
- Caution:** The shunt might not be preloaded correctly, if you notice an increased resistance already clearly before 1.5-2 cm to the marking or if you do not notice any resistance at all. DO NOT use and discard the product. If the force required to push the shunt in the cannula is inappropriate, the product must not be used.

*Figure 8: Insert the ejector through the cannula handle*

### C. Application of the cannula:

12. Select the puncture site and plan the route of the cannula to the fetus under ultrasound.
- Warning:** Make sure that you do not injure any organs (umbilical cord, etc.) with the cannula during the application.
- Warning:** Consider the length of the shunt (25 mm). When placing the shunt, ensure that the shunt ends (in the amniotic sac and fetal bladder) are exposed without damaging surrounding tissue / organs (e.g. umbilical cord) or impairing the drainage effect by kinking.
- Warning:** Consider the distance that needs to be bridged between the fetus abdominal wall and the fetus bladder in order to maintain a drainage effect. With growing fetus these dimensions might change.
13. Disinfect the area around the planned puncture site and, if necessary, apply local anesthetic in the area of the puncture site.
14. If necessary, make a stab incision of the puncture site with a scalpel to facilitate the penetration of the skin.
15. Remove protective tube (3) from the application cannula.
- Warning:** The cannula has a very sharp tip. There is a risk of injury.

16. Check the tip of the cannula if the shunt had been properly preloaded.


Tip is free: OK

Shunt is visible: Do not use and discard the product!

17. Puncture the fetal bladder under ultrasound (see *Figure 9*).

**Warning:** Be careful not to advance the ejector (7) during the puncture. Advancing the ejector may complicate or prevent the puncture of the fetal bladder.

*Figure 9: Application of the cannula in the fetal bladder*


#### D. Unfolding the shunt in the fetal bladder:

18. To unfold the distal brim of the shunt in the fetus, advance the ejector (7) to the second mark (b) (see *Figure 10*).

19. Check the position of the half-ejected shunt with ultrasound and, if necessary, correct the position of the shunt by moving the entire system (see *Figure 11*).

**Caution:** DO NOT move the cannula (2) and ejector (7) against each other. The mark (b) of the ejector remains in a constant position to the handle of the cannula (see *Figure 11*).

**Note:** The half-ejected shunt is approximately 1 - 1.5 cm outside of the cannula.


**Figure 10: Unfolding the distal brim of the shunt in the fetal bladder**

**Figure 11: Correcting the shunt position in the fetus**

#### **E. Unfolding the shunt in the amniotic sac:**

20. Carefully move the cannula handle (2) towards the proximal end of the ejector (7) until it meets the handle of the ejector (see *Figure 12*). Make sure that the ejector does not move during this step.

**Figure 12: Unfolding the proximal brim of the shunt in the amniotic sac**

The shunt is fully released and the second brim is unfolded if the cannula handle meets the handle of the ejector.

**Warning:** Before you fully release and place the shunt, reassure that the cannula tip is in the amniotic fluid (outside of the fetus and not in the uterine wall) and that the proximal end of the shunt will be securely placed in the amniotic fluid, not in the uterine wall. Consider shunt length of 25 mm.

21. Check the position and drainage performance of the shunt after placement under ultrasound.

*Image courtesy of Prof. Dr. med. Ulrich Gembruch, University Hospital Bonn*

#### **F. Disposal:**

22. Dispose of the product properly in a suitable container (follow internal regulations). Dress the wound.

#### **G. Follow-up:**

Close and frequent monitoring of the drainage performance under ultrasound (fetal bladder size in the normal range) must be carried out.

#### **H. Postnatal removal:**

Remove the shunt. The postnatal removal of the shunt can be associated with surgery.

**Warning:** Dislocated shunts should be surgically removed from the body of the newborn and / or the mother.


#### **Warning:**

The company SOMATEX does not assume any liability for the use of this product or its components in case of re-sterilization or reuse. This product may not be reused after a single application. The quality of the materials, coats and adhesive joints could degrade. Safe use is not guaranteed any longer. The product that is already used once is not designed for the required cleaning and sterilization processes. The sterility of the reprocessed disposable products is therefore not guaranteed. The risk of unwanted injuries and infections, especially cross-infections between patient and medical staff inappropriately increases.

#### **Storage Instructions:**

Store in a dry and cool area at a temperature of 5 – 30 °C (41 - 86 °F).  
Keep away from sunlight and heat.

**Any serious events that occur in relation to the product should be reported to SOMATEX Medical Technologies GmbH as well as the competent national authority.**

0482

**SOMATEX®**  
MEDICAL TECHNOLOGIES GMBH

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