# **Bone Marrow**


**REF**

180910

180935

180960

180965

**GEBRAUCHSANWEISUNG**  
**INSTRUCTIONS FOR USE**  
**MODE D'UTILISATION**  
**ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ**  
**ISTRUZIONI PER L'USO**  
**INSTRUÇÕES DE UTILIZAÇÃO**  
**MODO DE EMPLEO**  
**NÁVOD K POUŽITÍ**  
**UPUTE ZA UPORABU**  
**INŠTRUKCIE NA POUŽÍVANIE**

**INHALT / CONTENT**

Table des matières / Περιεχόμενα / Contenido / Índice/ Contenido / Obsah / Sadržaj:

**INHALT / CONTENT** ..... 2

**DEUTSCH** ..... 4

**ENGLISH**..... 6

**FRANÇAIS** ..... 8

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**ITALIANO**..... 12

**PORTUGUÊS**..... 14

**ESPAÑOL** ..... 16

**ČESKY** ..... 18

**HRVATSKI** ..... 20


**SLOVENČINA** ..... 22

**SYMBOLE / SYMBOLS**..... 24

**INFO** ..... 27

**ABBILDUNGEN / PICTURES**

Figures / Σχήματα / Figure / Imagens / Imágenes / Obrázky / Slike


Abbildung 1 / Figure 1 / Figure 1 / Σχήμα 1 / Figura 1 / Imagen 1 / Obrázek 1 / Slika 1


Abbildung 2 / Figure 2 / Figure 2 / Σχήμα 2 / Figura 2 / Imagen 2 / Obrázek 2 / Slika 2


Abbildung 3 / Figure 3 / Figure 3 / Σχήμα 3 / Figura 3 / Imagen 3 / Obrázek 3 / Slika 3


Abbildung 4 / Figure 4 / Figure 4 / Σχήμα 4 / Figura 4 / Imagen 4 / Obrázek 4 / Slika 4


Abbildung 5 / Figure 5 / Figure 5 / Σχήμα 5 / Figura 5 / Imagen 5 / Obrázek 5 / Slika 5

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*Keep until package is empty*

**Intended use:**

The **Bone Marrow Biopsy needle** is used for the puncture and histological extraction of tissue from the bone marrow in the iliac crest region.

**Contraindications:**

- Haemophilia or similar blood coagulation disorders
- Infection at the site of puncture

**Possible known complications:**

- Secondary haemorrhage and haematoma
- Infection
- Injury of the skin, soft tissues, nerves or adjacent organs
- Impairment of respiratory functions by the use of analgesics or sedatives

**Product description:**

The biopsy needle consists of a stylet with a fixed cap (1), a cannula with a grip (2) and a female Luer Lock connector for syringes, a depth stopper (3) to mark the depth of penetration, an insertion aid (4) and an ejector (5) to push out the tissue sample. The ejector has ring markings (6). The additional cap (7) helps to prevent further blood leakage.

The design follows the principles of the standard Jamshidi biopsy.

**Bone Marrow** is an advance in histological puncture methods, with an inner securing edge at the cannula tip. This feature supports specimen fixation inside the cannula.

**Warnings:**

- The products are **not** intended for **sternal** bone marrow aspiration.
- The products are not suitable for the injection of bone cement.
- For the safe and effective application of bone marrow puncture systems, the physician must have the relevant knowledge, experience and training in the use of this technique on patients.
- Follow aseptic procedures! Deep wound infections may occur as severe subsequent complications following the application, the removal of which may require surgical intervention.
- The products should be used in a clinical environment only.
- Sterile disposable gloves are required.
- Caution, risk of injury! The individual components of these products are very sharp and can easily cause injuries.
- These instructions are not intended as a description of interventional or surgical techniques. It is the responsibility of the attending physician to evaluate and to consider the possibility of the safe application of these products and techniques for each individual patient.
- Do not use if packaging is damaged or open.
- Never use the product after the sterility expiration date indicated on the packaging.
- The products are intended explicitly for single use only. The product may not be reused after a single application. The quality of the materials, coatings and adhesive bonds may deteriorate. A secure application is no longer guaranteed. The product is not designed for cleaning and sterilisation processes after a single use. The sterility of reprocessed single-use is therefore not guaranteed. The risk of unintentional injuries and infections, especially cross-infection in patients and medical staff considerably increases. The company SOMATEX® Medical Technologies holds no liability in the case of re-sterilisation and/or re-use of any device.
- The products are not designed for MRI safety areas. There is a risk of injury!
- It is possible, that no or too little of the biopsy material is extracted so that the application has to be repeated.

**Preparation:**

- The area surrounding the puncture site should be covered to create aseptic conditions.
- Appropriate measures to secure the biopsy sample (biopsy material slide, sample container) for haemopathological evaluation should be taken prior to the intervention.

**Biopsy application:*****(Follow the illustrations!)***

1. Perform subcutaneous and periosteal local anaesthesia at the puncture site.
2. If necessary, use a sharp scalpel to make an incision at the puncture site to facilitate the penetration of the skin.
3. Remove the protective tube from the cannula. The depth stopper (3) should remain on the cannula.
4. **Note:** Prior to the puncture, the insertion depth stopper (3) can be adjusted on the cannula shaft for visual control of the penetration depth (see figure 1).
5. Position the puncture needle (2) and push it forward until you feel a noticeable resistance indicating that you reached the bone tissue.
6. Penetrate the needle into the iliac crest in a targeted manner using continuous pressure while rotating (right/left rotation).
7. Upon reaching the medullary canal, the resistance will decrease as it penetrates. Then, pull the trocar out of the puncture needle (see figure 2) by lifting off the cover (1).
8. Optionally, blood can now be aspirated from the medullary canal using a syringe.
9. To achieve better handling and to avoid the further leakage of blood from the Luer-Lock connector, place the additional cap (7) on the needle handle (see figure 2).
10. Now push the needle further into the medullar canal by means of rotation movements until the required biopsy amount is obtained (see figure 3).  
**Note:** The length of the biopsy sample can be checked by inserting the ejector (5). For this purpose, remove the additional cap (7) from the handle and insert the ejector (5) carefully into the cannula until you feel a noticeable resistance. The cm markings (6) indicate the length of the biopsy sample within the puncture needle (see figure 4).
11. Move the needle with a suitable strain back and forth to separate the sample from the marrow.  
**Caution: Avoid unsuitable, excessive force to prevent the possible bending or breaking of the needle.**
12. Then remove the needle from the bone.
13. Treat the wound.

**Collection of the biopsy sample from the needle after the puncture:**

1. Remove the replacement cap (7) from the puncture needle (2).
2. Place the insertion aid (4) on the tip of the puncture needle in the direction of the arrow.
3. Insert the ejector (5) into the conical opening of the insertion aid and carefully remove the biopsy sample from the cannula in a distal to proximal direction (see figure 5).
4. Dispose of the used puncture needle in the appropriate manner.

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SOMATEX Medical Technologies GmbH