
Lung Marker System

Lungenmarkersystem

REF 272012 272015

GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MODE D'UTILISATION
ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ
ISTRUZIONI PER L'USO
INSTRUÇÕES DE UTILIZAÇÃO
MODO DE EMPLEO
NÁVOD K POUŽITÍ
NAUDOJIMO INSTRUKCIJOS
BRUGSANVISNING
KÄYTTÖOHJEET
GEBRUIKSAANWIJZING
INSTRUKSJONER FOR BRUK
INSTRUKCJA UŻYCIA
ANVÄNDNINGSSINSTRUKTIONER

INHALT / CONTENT

INHALT / CONTENT	2
DEUTSCH	3
ENGLISH.....	5
FRANÇAIS	7
ΕΛΛΗΝΙΚΑ	9
ITALIANO.....	11
PORTUGUÊS.....	13
ESPAÑOL	15
ČESKY	17
HRVATSKI	19
LIETUVIŲ	21
DANISH	23
SUOMI	25
NEDERLANDS.....	27
NORSK.....	29
POLSKI	31
SVENSKA	33
SYMBOLE / SYMBOLS	35
INFO	40

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
ENGLISH

**Intended use:**

The SOMATEX® Lung Marker System is used for the pre-operative marking of intrapulmonary, non-palpable pulmonary nodules. By using this system, localization of small lesions is facilitated by means of a minimally invasive, endoscopic technique.

Note:

This product was expressly manufactured for single use and may not be reused after its initial use. SOMATEX® declines any responsibility in the case of product resterilization or reuse. The product is not designed for the essential cleaning and sterilization processes. Sterility of the reprocessed single-use product is therefore not guaranteed. In addition, the quality of the material, coating and adhesive bonds may be impaired, with safe use no longer be guaranteed. The risk of unintended injuries and infection, particularly cross-infection in patients and medical personnel, increases disproportionately in this case.

This product may only be used by qualified medical personnel.

Contraindications:

Absolute contraindications: increased bleeding propensity, no suitable puncture site localizable

Relative contraindications: intake of anticoagulants, insufficient lung reserve, cardiac insufficiency, pulmonary hypertension, pronounced emphysema, chronic respiratory insufficiency, severe dyspnoea, severe pleural adhesions

Possible known complications: slight bleeding at the puncture site and needle path, pneumothorax, wire dislocation, pain, emphysema, dyspnoea, haemoptysis

Size of the wire coil: The diameter of the spiral wire (4) is (8 ± 1) mm with a height of (6 ± 1) mm.

The following steps must be adhered prior to use:

1. Ensure that the packaging is unopened and/or undamaged and check the expiration date. In the case of damaged packaging, the product may not be used.
2. Open the packaging
Note: On delivery, the spiral part of the wire (4) is outside, and the straight wire shaft inside the loading cannula (3). The fixation part (2) is premounted on the free end of the wire.
3. Prior to the application, the spiral wire must be completely retracted into the loading cannula by pulling on the fixation part (2).
4. The fixation part (2) can be moved along the wire or removed after loosening the screw cap. The system is thus prepared for its application. The intended placement of the wire should be agreed with the thorax surgeon beforehand.

Note: The wire spiral should preferably be placed behind the lesion or at its periphery.

Description of use:

5. Local disinfection and administration of local anaesthetic at the puncture site
6. Use a scalpel to carry out an incision at the puncture site to facilitate penetration of the skin, if necessary.
7. Insert puncture needle (1) - guide the puncture needle to the target area under imaging guidance.
Note: The position and orientation of the spiral should be documented and transmitted to the thorax surgeon.
8. After reaching the target area, only the mandrel of the puncture needle (1) should be removed.
9. Completely insert the loading cannula with the retracted wire into the positioned puncture needle.
10. Carefully insert the wire from proximal to distal in the system. This will result in the deployment and tissue fixation of the spiral in the target area. Hold the wire approx. 1 cm from the cannula hub and slide it gradually into the cannula. The complete deployment and correct positioning of the wire marker is to be verified by selected imaging. If a position correction is necessary, the wire can be retracted in the chest cannula and repositioned, as described in step 3.
11. After placement of the wire, remove the fixation part (2) from the wire shaft by loosening the screw cap.
12. The needle system can now be removed through the wire, which remains in the patient.
13. The fixation part can now be deployed again with the flat side facing the surface of the skin. It can then be fixed by tightening the screw cap.
Note: When mounting of the fixation part, the respiratory movement of the chest should be considered to reduce the risk of wire dislocation.
14. Treat the wound. The outer part of the wire during is to be properly secured and protected transthoracically (tape, dressing).
Note: To reduce the possibility of unwanted complications (pneumothorax, wire dislocation) the period of time up to resection of the lesion should be minimized.
15. All components of the localization system should be appropriately disposed in a designated container.

Note to the surgeon:

The "hanging" lesion can be carefully pulled outward intraoperatively with controlled force, verified endoscopically and resected from the tissue with appropriate instruments. Use the specimen pouch to remove the resected tissue, along with the intrathoracic portion of the marker wire, through the appropriate trocar puncture site. Where appropriate, the outer part of the wire can be severed extrathoracically beforehand.

Note:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is resident.

Directions for storage:


Keep away from sunlight and heat. Keep dry.

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SOMATEX Medical Technologies GmbH