
# Q-Wire System

REF

271462 271463 271464 271465

**GEBRAUCHSANWEISUNG**  
INSTRUCTIONS FOR USE  
MODE D'UTILISATION  
ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ  
ISTRUZIONI PER L'USO  
INSTRUÇÕES DE UTILIZAÇÃO  
MODO DE EMPLEO  
INSTRUCȚIUNI DE UTILIZARE  
NÁVOD K POUŽITÍ  
UPUTE ZA UPORABU  
HASZNÁLATI ÚTMUTATÓ  
GEBRUIKSAANWIJZING  
INSTRUKSJONER FOR BRUK  
INSTRUKCJA UŻYCIA  
INŠTRUKCIE NA POUŽÍVANIE  
ANVÄNDNINGSSINSTRUKTIONER  
BRUGSANVISNING

**INHALT / CONTENT**

CONTENU / ΠΕΡΙΕΧΟΜΕΝΟ / SODDISFARE / CONTENTE / CONTENIDO / OBSAH / SADRŽAJ / TARTALOM /  
 INHOUD / INNHOLD / ZAWARTOŚĆ / INNEHÅLL / INDHOLD:

<b>INHALT / CONTENT</b> .....	<b>2</b>
<b>DEUTSCH</b> .....	<b>3</b>
<b>ENGLISH</b> .....	<b>5</b>
<b>FRANÇAIS</b> .....	<b>7</b>
<b>ΕΛΛΗΝΙΚΑ</b> .....	<b>9</b>
<b>ITALIANO</b> .....	<b>11</b>
<b>PORTUGUÊS</b> .....	<b>13</b>
<b>ROMÂNĂ</b> .....	<b>15</b>
<b>ESPAÑOL</b> .....	<b>17</b>
<b>ČESTINA</b> .....	<b>19</b>
<b>HRVATSKI</b> .....	<b>21</b>
<b>MAGYAR</b> .....	<b>23</b>
<b>NEDERLANDS</b> .....	<b>25</b>
<b>NORSK</b> .....	<b>27</b>
<b>POLSKI</b> .....	<b>29</b>
<b>SLOVENČINA</b> .....	<b>31</b>
<b>SVENSKA</b> .....	<b>33</b>
<b>DANSK</b> .....	<b>35</b>
<b>SYMBOLE / SYMBOLS</b> .....	<b>37</b>
<b>INFO</b> .....	<b>42</b>

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## Instructions for use

**Important Information:**

Read this instruction manual thoroughly and be familiar with its contents prior to use. Failure to read the entire manual and familiarize yourself with all instructions before using the **Q-Wire System** is unsafe and can result in life threatening or severe injury to the patient or user and to damage or malfunction of the device.

**Indications:**

The product is intended as a preoperative marker of non-palpable suspected breast lesions to facilitate the intraoperative localization of the findings by the surgeon.

It is not indicated to be used with magnetic resonance imaging (MRI) techniques.

**Contraindications:**

All contra-indications applicable to the relevant area of application, as known according to the rules of the art of medicine and anticipated for the use of cannulas and marker systems for preoperative marking of breast lesions, shall apply.

**Possible known complications:**


Wire dislocation, accidental cutting of wire, wire break, bleedings, infections, aesthetic complications

**Warnings:**

- Only qualified physicians with knowledge, experience and training in percutaneous soft tissue marking should use the **Q-Wire System**.
- This manual does not include descriptions or instructions for surgical techniques. It is the responsibility of the physician performing any procedure to determine the appropriateness of the type of procedure to be performed and the use of this product and to determine the specific technique for each patient.
- DO NOT use the system in patients with breast implants.
- The **Q-Wire System** tissue site marking system is only sterile if used before the expiration date and if package is unopened and undamaged. DO NOT use after the expiration date or if package is open or damaged.
- Single patient use only. DO NOT reuse or resterilize.
- The company SOMATEX does not assume any liability for the use of this product or its components in case of re-sterilization or reuse. This product may not be reused after a single application. The quality of the materials, coats and adhesive joints could degrade. Safe use is not guaranteed any longer. The product that is already used once is not designed for the required cleaning and sterilization processes. The sterility of the reprocessed disposable products is therefore not guaranteed. The risk of unwanted injuries and infections, especially cross-infections between patient and medical staff inappropriately increases.
- The product is NOT suitable for MRI use (magnetic resonance imaging)! Danger of injury!
- The product is NOT suitable in direct contact with electrosurgical instruments! Danger of spark discharging and injury!
- Care must be taken, when using electrosurgical instruments: the product can be damaged in direct contact. If damage to the marking wire occurs, a follow-up check and, if necessary, removal of wire fragments remaining in the body is required.

**Precautions:**


- Cannula tip is sharp. Use care especially while unpacking the cannula.
- Make sure that the marker wire is located inside of the cannula and that the curved part of the wire do not protrude out of the tip of the cannula. Both circular marks (1) of the wire must be visible. (see Picture 1)
- **Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician.**


Picture 1

**Device Description:**

The **Q-Wire System** is a sterile, single use, tissue site marking system consisting of a nickel-titanium marker wire (3) and a puncture needle (2) and 1 cm depth marks to aid in cannula placement (see Picture 2). The puncture needle is designed with a bevelled tip for convenient introduction. The marker wire is designed with two different markings, the first (1) marker allows the use to confirm that the curved part of the wire do not protrude out of the tip of the cannula prior to deployment. To avoid the risk of a change of the wire position in the distal direction a wire fixation (4) is usable.

Picture 2

Picture 3

**Directions for Use:**

1. Prior to opening the package, make sure that the package has not been opened and/or damaged. In addition, check the sterilization expiration date.
2. Open the package.
3. Disinfect the puncture area and administer local anesthetic in the insertion area.
4. Use the scalpel to carry out an incision in the puncture site to facilitate the insertion into the skin.
5. Before starting the puncture make sure that the curved part of the wire do not protrude out of the tip of the cannula. Both circular marks (1) of the wire must be visible.
6. Needle insertion: insert the needle under sonographic or mammographic guidance into the mamma such that the tip of the needle reaches the tumor itself or the tumor region.
7. After positioning the needle in the target area, the marker wire can be placed in the tumor area by holding the needle position and pushing the wire forward into the cannula. The marking wire tip is fully extended and in a curved shape once the second mark (1) on the marking wire gets invisible. (Picture 3)

**Note:** The wire should be pushed into the cannula only until the last visible marking on the wire (1) has been inserted into the cannula attachment!

**IMPORTANT:** Do not advance wire into the cannula any further past the second marking on the wire!

8. Prior to the removal of the needle, a proper positioning of the rolled wire can be inspected. Should the position be other than optimal, the wire can be withdrawn into the cannula. After rectifying the needle position, the wire can be deployed again according to the description in points 5 and 6 for preoperative marking.
9. Carefully remove the needle (2) from the mamma. If possible, bring the patient in an upright position. Relax the compression of the breast, if necessary. Slide the opened fixation part over the wire almost to the surface of the skin of the breast (**Caution:** The breast must be decompressed when this is done. There must still be some space between skin and fixation part). The wire must lie between the brackets (4a) of the fixation part. The wire is fixed through the closure being locked in place (4b). This measure prevents the risk of change in position of the wire along the distal end.
10. Treat the wound.
11. Dispose the product properly in a suitable container.

**Storage Instructions:**

Store in a dark, dry place with a temperature between 5° – 30°C (41° - 86° F).

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**SOMATEX Medical Technologies GmbH**