



# Biopsie-Handy

## Halbautomatisches Biopsiesystem

### Vorteile im Überblick:

- Bequeme Einhandbedienung
- Geringes Gewicht
- Ultraschall- und CT-gestützt einsetzbar
- Atraumatische Nadelspitze
- cm-Markierung
- optimal mit koaxialer Schleuse (koaxiale Punktionskanüle) erhältlich
- Auch als MR-kompatibles System verfügbar



Biopsie-Handy			
Artikelnr.	Gauge	Durchmesser	Länge
900 118	20	0,90 mm	100 mm
900 120	20	0,90 mm	150 mm
900 122	20	0,90 mm	200 mm
900 126	18	1,20 mm	100 mm
900 128	18	1,20 mm	150 mm
900 130	18	1,20 mm	200 mm
900 134	16	1,60 mm	100 mm
900 136	16	1,60 mm	150 mm
900 138	16	1,60 mm	200 mm
900 142	14	2,10 mm	100 mm
900 144	14	2,10 mm	150 mm
900 146	14	2,10 mm	200 mm

Passende Punktionschleuse



Koaxiale Punktionskanüle (Punktionschleuse)			
Artikelnr.	Gauge	Durchmesser	Länge
170 251	18	1,20 mm	50 mm
170 253	18	1,20 mm	100 mm
170 255	18	1,20 mm	150 mm
170 257	17	1,50 mm	50 mm
170 259	17	1,50 mm	100 mm
170 261	17	1,50 mm	150 mm
170 263	15	2,00 mm	50 mm
170 265	15	2,00 mm	100 mm
170 267	15	2,00 mm	150 mm
170 269	13	2,40 mm	50 mm
170 271	13	2,40 mm	100 mm
170 273	13	2,40 mm	150 mm

# Biopsy-Handy

## Semi-Automatic Biopsy Device

### Advantages at a glance:

- Easy one-hand operation
- Light weight
- Ultrasound and CT-compatible
- Atraumatic needle tip
- cm-markings
- Optionally available with coaxial sheath (coaxial puncture needle)
- For application with MRI control the MRI-system is available



Biopsy Handy			
Order No.	Gauge	Diameter	Length
900 118	20	0,90 mm	100 mm
900 120	20	0,90 mm	150 mm
900 122	20	0,90 mm	200 mm
900 126	18	1,20 mm	100 mm
900 128	18	1,20 mm	150 mm
900 130	18	1,20 mm	200 mm
900 134	16	1,60 mm	100 mm
900 136	16	1,60 mm	150 mm
900 138	16	1,60 mm	200 mm
900 142	14	2,10 mm	100 mm
900 144	14	2,10 mm	150 mm
900 146	14	2,10 mm	200 mm

Fitting puncture Sheath

Coaxial Puncture Needle (Biopsy Sheath)			
Order No.	Gauge	Diameter	Length
170 251	18	1,20 mm	50 mm
170 253	18	1,20 mm	100 mm
170 255	18	1,20 mm	150 mm
170 257	17	1,50 mm	50 mm
170 259	17	1,50 mm	100 mm
170 261	17	1,50 mm	150 mm
170 263	15	2,00 mm	50 mm
170 265	15	2,00 mm	100 mm
170 267	15	2,00 mm	150 mm
170 269	13	2,40 mm	50 mm
170 271	13	2,40 mm	100 mm
170 273	13	2,40 mm	150 mm