


# SurgLite

## Surgical illumination made easy

SOMATEX® SurgLite provides excellent illumination of the surgical field without any time- and cost-consuming preparation. The flexible material as well as different sizes and run-times ensure maximum reliability. SurgLite can be easily fitted and attached to a variety of retractors before or during surgery.

*Possible applications:*

*Breast surgery (e.g. subcutaneous mastectomy),  
general surgery, plastic surgery,  
orthopedic surgery, cardiothoracic surgery,  
urology, ENT (thyroid gland)*



- **Time and cost savings** thanks to flexible application and no need of sterilization, repair or maintenance
- **Optimal visual support** even in cavities that are difficult to access
- **Comfortable handling** due to easy and quick fixation, flexible material, no light guides to disturb the surgeon, no heat generation due to high efficiency LEDs

Article No.	SL0200-15	SL0100-15	SL0200-40	SL0100-40	SL0200-60	SL0100-60
Strip length	200 mm	100 mm	200 mm	100 mm	200 mm	100 mm
Optimal run-time	15 min	15 min	40 min	40 min	60 min	60 min
Packaging unit: 10 pieces						

SOMATEX Medical Technologies GmbH  
Hohenzollerndamm 150 / 151 · 14199 Berlin · Germany  
Tel.: +49 (0) 30 31982 25-00 · Fax: +49 (0) 30 31982 25-99  
info@somatex.com · www.somatex.com

**SOMATEX**  
MEDICAL TECHNOLOGIES